

Process Analytics

Measuring Solutions



The World of Process Analytics

Learn more on our Website

Knowledge Base
Latest Innovations & Software Updates
Manuals & Specifications
Application Notes
Quality & Regulatory Certificates



www.hamiltoncompany.com/process-analytics

Keep yourself updated


Follow us on the LinkedIn Channel



www.linkedin.com/showcase/hamilton-process-analytics

Process Analytics Measuring Solutions





See more on page 4



Single Use Sensors

See more on page 5



CO₂NTROL

See more on page 6



VisiFerm RS485

See more on page 6




Table of Contents

Highlights4
Innovations 8
Arc Family 12
System Installations 14
Sensors 16
pH16
ORP
Conductivity 54
Cell Density 66
CO ₂
DO 78
DuraCal pH Buffers 100
Conductivity Standards 102
Oxygen Accessories 104
Electrolytes and Solutions 105
Connectivity106
Cables 108
Arc Accessories 117
Hamilton Customized Products 119
Transmitter 120
H100 120
H220X 124
Housings 126
Sensor Comparison 160
Safety First162
Alphabetical Index 164

Highlights

CELL DENSITY

On-Line Data Real-Time Decisions





Permittivity measurements are the most reliable method of monitoring Viable Cell Density (VCD). This measurement is immediately affected by changes in Viable Cell Density and can be used to time process-specific actions for maximum yield. Permittivity can also be used to detect changes in cell physiology and is the most immediate method for determining the beginning of the cell death phase





Next Generation Total Cell Density Measurement with Dencytee Arc

With the Dencytee Arc sensor, Hamilton now offers a new generation of in-line total cell density biomass measurement. We have taken the measurement technology to the next level and combined the advantages of transmission and reflection measurement. By upgrading to two detectors, higher measurement resolution can be achieved. This results in higher reliability that can be used in both low and high cell concentrations.





SINGLE USE One Vendor All Measurements

Hamilton has worked closely with single-use (SU) equipment manufacturers to understand the market needs in order to adapt measurement technologies from reusable sensors because all applications have their own requirements. The Hamilton SU sensors offer the known high accuracy of traditional sensors even after gamma irradiation and dry storage. The SU portfolio offers sensing elements as well as a wide variety of possible connections to transmitters and controllers. Arc modules are also available for easy integration of digital signals and allow, in combination with the ArcAir app, to benefit from the Arc technology. Thus calibration data provided on a label can easily be scanned and the sensors are ready to be used with seconds.



VisiFerm DO SU Family

Reliable Dissolved Oxygen Measurement

The Hamilton VisiFerm DO SU sensor systems are available in a wide application range for bag and rigid containers. Various mechanical connections in the vessel are available with a single-use sensor element and reusable electronic for a cost effective application. The new single use optical dissolved oxygen sensor offers a reliable and comparable measurement to existing reuseable probes.



OneFerm pH Family

High Performance pH Measurement

The Hamilton OneFerm pH sensor is a single use glass electrode in order to ensure a wide measuring range, and a very low drift, even after dry storage and wet-in time. Sensors are available in various lengths and electrical connections so that the pH measurement can benefit from the Arc technology.



Incyte SU Family

Monitoring Viable Cell Density

Online cell density measurement is essential to ensure reliable processes, especially for long running, i.e. perfusion. Online data provides continuous information in order to optimize control and yield.

Conducell SU Family

Conductivity Measurement In Bags

The Conducell SU Family allows measurements in a wide conductivity range in SU applications.



HIGHLIGHTS

Get Co₂ntrol

Solid-State Optical DCO₂ Sensor

Though DCO_2 is commonly recognised as a critical process parameter in biopharmaceutical production, the measurement technology has not really changed a lot. In fact, all in-line sensors available on the market until now are based on the indirect Severinghaus measuring principle – a technology that is more than 50 years old and prone to measurement errors and high maintenance.

It was clear that Hamilton Process Analytics would take on the challenge to develop a new type of sensor that would combine real-time control together with reliability and cost efficiency.

We are now more than proud to present you CO₂NTROL – our new solid state sensor that directly measures DCO₂ and provides maintenance free, real-time, and in-line control of this critical process parameter.

Find all details about our new sensor on pages 74 to 77.



Intuitive Sensor Management

The ArcAir App: One Tool for Sensor Management & Ready for GMP Compliance



- Wireless configuration and calibration
- Common interface for mobile, tablet, and PC
- Automated validation and documentation
- ▶ Ready for compliance with FDA CFR 21 Part 11 and Eudralex Volume 4 Annex 11













VisiFerm RS485

The Sensor for Biopharma

Hamilton fully redesigned the sensor electronics and optical cap to create the most robust VisiFerm dissolved oxygen sensor yet. Upgrading both key components allowed the VisiFerm RS485 to have less frequent need for calibration, less measurement error, and longer lifetime than previous optical sensor thechnologies.

- ▶ 80% Fewer Calibrations
- 3x Longer Cap Life
- ► 50% Longer Sensor Life
- Supply security for the upcoming decade

VisiTrace RS485

Trace Level DO Measurement

Hamilton fully redesigned the sensor electronics to create the most robust VisiTrace sensor yet. VisiTrace RS485 is designed for trace measurement from 1 to 2,000 ppb and stable against active chlorine and chlorine dioxide.



Arc Modbus OPC Converter

Easy Integration in SCADA

The Arc Modbus OPC converter is designed to connect Arc Modbus to OPC UA. It is perfect to be used in R&D areas for the integration via Ethernet into SCADA systems and supports all VP8 Arc sensors. Up to 4 sensors in parallel can be used. The conversion script is pre-installed and the converter ready-to use.



INNOVATIONS

Beyond Process Analytics

Hamilton's electrochemical and optical sensors are the solution for process analytical measurement systems, characterized by proven quality and outstanding performance. Offering measurement parameter solutions in pH, ORP, dissolved oxygen and conductivity, our sensors and accessories are backed by over 50 years of engineering and manufacturing expertise in innovative design.



pH Glasses

Measurement Accuracy in Various Applications

Measurement stability and sensor lifetime in various environments requires different pH glasses.

Our high performance glasses, the PHI and the HB glass, were developed to withstand frequent steam sterilization, autoclaving and CIP cleaning using hot caustics. PHI and HB glass provide the lowest drift and show almost no shift after sterilization and cleaning procedures.

The H glass has excellent aging characteristics and offers stable readings even in samples with low water content such as anhydrous or only partially aqueous solutions. The low alkali error of H glass means accurate measurements even at high pH or high operating temperatures. HF glass ensures the longest possible lifetime in low temperature processes and processes containing hydrofluoric acid.

Foodlyte

Biocompatible Reference Electrolyte

The Foodlyte electrolyte was specifically developed for the needs of the biotechnology, pharmaceutical and food industries. It's based on food ingredients and the perfect electrolyte for applications where non-toxicity is mandatory. Foodlyte is taste-, odor- and harmless for microorganisms.

The biocompatibility is approved by MDT¹ according to EN ISO 10993-5² and USP 31, 2008 Chapter 87³ and according to international GLP⁴ guidelines.



- 1 Medical Device Testing GmbH Ochsenhauser
- 2 Biological evaluation of medical devices -- Part 5: Tests for in vitro cytotoxicity
- 3 Biological Activity Tests, In Vitro
- 4 Good Laboratory Practice



Single Pore Concept

The never-clog Liquid Junction

A Single Pore is an open liquid junction and an alternative to diaphragms. Instead of many tiny pores in a ceramic diaphragm, a single pore, about 2000 times larger in diameter, is used. This concept provides a direct contact between reference electrode and sample. In combination with the bigger diameter this liquid junction can hardly be clogged. The Single Pore results in a faster response time, more accurate readings and prevents reference poisoning.

Note: The PTB (Physikalisch-Technische Bundesanstalt = Physical Technical Federal Institute) in Braunschweig, Germany, determined the Single Pore pH electrode to be the most accurate laboratory electrode. Further information can be found in "Traceability of pH measurement" by Petra Spitzer; ISBN 3-89429-877-4 or ISSN 0947-7063.

Polisolve Plus

Most innovative Polymer Reference Electrolyte

Hamilton has designed innovative Polisolve Plus polymer electrolyte sensors that cover the full pH range, a wide temperature range and withstand reference poisoning for an extended lifetime. It's also stable against most organic solvents and free of toxic acrylamide.

When Polisolve Plus and Single Pore concepts are combined the result is a Polilyte Plus sensor for a wide range of applications as well as a problem solver for difficult applications.

- Industrial waste water
- Hot sugar juice
- Samples containing color pigments
- Oily samples

The combination leads to more stable reference signals and minimized diffusion potentials. Polisolve Plus represents a significant contribution to long lasting pH sensors.



Conductivity Standards

Certified and Traceable

Hamilton was the first to offer conductivity standards at 1.3 and 5 μ S/cm with a certified accuracy of ±1% and a durability of 1.5 or 3 years. All conductivity standards exhibit a previously unknown level of stability which has been confirmed by measurements done by the PTB¹. Governmental metrological institutes that deal with measurement of electrolytic conductivity have become aware of these standards, and the composition of these standards is patented. The measurement procedure for determining conductivity has been developed in collaboration with the DFM². Each batch is certified by the DFM. In an inter-laboratory test among prestigious European metrological institutes (PTB, DFM, DAkkS³), Hamilton standards were used as a measurement solution.

1 PTB: Physikalisch-Technische Bundesanstalt, Braunschweig, Germany 2 DFM: Danish Institue of Fundamental Metrology, Lyngby, Denmark 3 DAkkS: Deutsche Akkreditierungsstelle, Wolfen, Germany





DuraCal pH Buffers

Easy Calibration with 5-Year Shelf Life

DuraCal pH buffers consist of a complete range of patented stable pH buffer solutions from pH 1.09 to pH 12.00. Hamilton guarantees that they will last for five years from the date of manufacture. The pH 9.21 and pH 10.01 buffers are even stable in air. High buffer capacities enable quick and stable calibrations.

Closed-loop traceability: In contrast with other manufacturers Hamilton has developed a "closed-loop" traceability. For users of DuraCal pH buffer solutions this means a unique level of reliability.

Top-down traceability: With Hamilton the pH value of the DuraCal buffer is determined by a comparison with two secondary reference solutions.

Bottom-up traceability: From each lot manufactured, a representative quantity is measured at DAkkS (Deutsche Akkreditierungsstelle, Wolfen, Germany). This ensures an external independent verification by an accredited institute. The DAkkS issues an official calibration certificate for every DuraCal batch manufactured.



VisiFerm DO

The most reliable Optical Dissolved Oxygen sensor in the Industry

The VisiFerm DO is the first optical dissolved oxygen (DO) process sensor for demanding applications in the pharmaceutical, biotechnology and beverage industries. The measuring principle is based on oxygen dependent quenching of the emitting light of a luminophore. Easy and fast to maintain, the multiple time-constraints caused by the use of electrochemical type DO sensors is eliminated. Decreased cost of ownership is further improved with an integrated sensor lifetime check that indicates when the sensor is in need of maintenance. A simple, replaceable cap rebuilds the sensor in seconds.

The optical measurement is independent from the flow and insensitive to CO₂. A special window behind the luminophore enables the sensor to withstand pressure hammers and spikes. Due to this design, the VisiFerm DO is suitable for inline measurement of dissolved oxygen in various processes.



The True Power

Intelligence Integrated

Hamilton Arc revolutionizes the integration of sensors by rethinking communication between sensors, end users and process control systems (PCS). The functionality of a traditional transmitter has been replaced by a microprocessor within the sensors head. Arc sensors communicate directly with the PCS through 4-20 mA standard and digital signals.

Arc sensors offer a fully compensated, converted digital and 4-20 mA signal directly to the process control system.

Fully compensated signal

- ► Temperature compensated
- ► E.g. Pressure, Salinity

Signal output

- Digital Modbus
- ► 4-20 mA analog
- ▶ Different parameter units (e.g. mV, ppb, %sat...)






Arc Intelligence

Wireless Communication & Calibration

Arc sensors provide full online wireless option for monitoring, configuration and calibration.



Laboratory Calibration




Complete Arc Sensor Portfolio




Analog Systems


Standard Measuring Loop



Measuring Loop in Hazardous Area




Measuring Loop in Pipe





Arc Systems

Skid System




Arc in R&D





The pH value is measured in most processes using a glass electrode. This pH glass forms a thin gel layer in aqueous solutions that is highly selective to H⁺ ions. The pH dependent potential of the gel layer is measured against a built-in reference electrode with a constant potential. This reference electrode may be a silver wire in contact with solid silver chloride or a calomel electrode.

In general, the pH value is a measure of the acidity or the basicity of an aqueous solution. In technical terms, pH is the negative logarithm of the activity of the solvated protons H⁺. It's mostly explained as the measure of the proton concentration which is correct for dilute aqueous solutions.



Polilyte Plus family



For more specifications see www.hamiltoncompany.com

Specifications

Measuring range
Process temperature

Pressure range (relative to ambient)

Hygienic aspects

pH glass

Electrolyte

Diaphragm

O-ring

Reference system

0 to 14 pH

CIP: HB, PHI SIP: H, HB, PHI

Polisolve Plus

Single Pore

EPDM: HB, PHI FKM: H, HF

Everef-L

See table on page 160/161 See table on page 160/161

Autoclavable: H, HB, PHI

See table on page 18

The outstanding success of the Polilyte Plus in chemical and wastewater applications gave the inspiration for transferring the good features to a whole family of sensors. The expanded portfolio widens the range of applications that can be covered.

All members have the same reference electrolyte Polisolve Plus, use the Single Pore technology but will have different pH glasses.

Benefits

- ► More applications with HB pH glass
- ► Better overview of the portfolio
- ► There's always at least one family member that suits the different applications
- ► Resistant against solvents, strong acids and bases

Typical applications

- Sugar industry
- Microelectronics
- ► Industrial wastewater
- Downstream processes
- Fermentation







How to choose the sensor	New sensor	pH glass	Electrolyte	Predecessor
HF in the media, low temperature	Polilyte Plus HF	HF	Polisolve Plus	ClaryTrode
Low conductivity	Polilyte Plus H	Н	Polisolve Plus	Polilyte HT
CIP, SIP, autoclavations, chemical robustness	Polilyte Plus PHI	PHI	Polisolve Plus	Polyclave
CIP, SIP, autoclavations, fast response time	Polilyte Plus HB	HB	Polisolve Plus	
High pressure	Polilyte Plus XP	Н	Polisolve Plus	Polilyte Plus XP

Ordering Information

242428	Basic n	umber =	Polilyte P	lus VP 12	0 (old Ref)
	Code	pH glas			
	1	Н			
	2	HB (not	for MS)		
	3	HF			
	4	PHI			
		Code	Electri	cal Conn	ector
		1	VP 🚱		
		2	S8 🚱		
		3	Arc		
		4	Memos	ens 😉	
			Code	a-leng	th (mm)
			1	120	
			2	225	
			3	325	
			4	360 (no	ot for Arc, MS only with H glass)
			5	425	
				Code	Temperature sensor
				1	Pt100 (VP) (not applicable for Arc)
				2	Pt1000 (VP) (not applicable for Arc)
	+	+	+	3	none (S8) or given (Memosens, Arc)
242428 -					← Order Code



Accessories

pH buffers see page **▶** 100

Cables see page ▶ 108

Housings see page **▶** 126

SENSORS **SENSORS**

EasyFerm Plus family



The EasyFerm Plus family of pH sensors is designed to withstand demanding applications in the Pharmaceutical and Chemical industries. All family members have the same reference electrolyte Phermlyte, the same type of diaphragm HP Coatramic but different pH glasses. The standard EasyFerm Plus, with its PHI glass, is directed at the BioPharm and Pharmaceutical industries because the glass has an excellent chemical robustness and provides best results in applications where sterilization either in an autoclave or an SIP is performed frequently. The new versions with the HB glass show a very fast recovery after CIP and SIP cycles leading to a shortened set-up time.

The LEVP (LE = Liquid Earth) versions have a stabilized sensor signal and an extended sensor diagnosis.



Did you know... that with a pre-pressurized reference system the life time of a sensor is extended?

Benefits

- ► Pre-pressurized reference electrolyte ensures a clog-free diaphragm
- ► Almost drift-free measurement
- ► Stable measurement signals after steam sterilization, autoclavation and CIP cleanings

Typical applications

- ► Industrial processes








How to choose the sensor	New sensor	pH glass	Electrolyte	Predecessor
CIP, SIP, autoclavations, chemical robustness	EasyFerm Plus PHI	PHI	Phermlyte	EasyFerm Plus
CIP, SIP, autoclavations, fast response time	EasyFerm Plus HB	НВ	Phermlyte	





Ordering Information

238633					
	Code	pH glas	SS		
	1		ommended	d pH glass	type)
	2	НВ			
		Code	Electric	al Conne	ctor
		1	VP 😉		
		2	S8 😥		
		3	Arc		
		4	Memose	ens 😉	
		5	K8 😉		
		6	LEVP (o	nly for 120	and 225 mm length) 69
			Code	a-lengt	h (mm)
			1	120	
			2	160	
			3	200	
			4	225	
			5	325	
			6	360 (not	for Arc and only PHI glass)
			7	425	
			8	275	
				Code	Temperature sensor
				1	Pt100 (VP, LEVP) (not applicable for Arc)
				2	Pt1000 (VP, LEVP) (not applicable for Arc)
	+	+	+	3	none (S8, K8) or given (Memosens, Arc)
238633 –					← Order Code



Accessories

pH buffers see page **▶** 100

Cables see page ▶ 108

Housings see page ▶ 126

EasyFerm Bio [family]



Specifications Measuring range 0 to 14 pH Process temperature 0 to 140 °C (Arc: analog 0 to 110 °C, digital 0 to 140 °C) 0 to 6 bar Pressure range (relative to ambient) Hygienic aspects Autoclavable, SIP, CIP HB, PHI pH glass Electrolyte Foodlyte Reference system Everef-F Diaphragm **HP** Coatramic O-ring Silicone

For more specifications see www.hamiltoncompany.con

The EasyFerm Bio family of pH sensors is designed for applications in the Pharmaceutical, Biotechnology and Food & Beverage industries. All family members have the same reference electrolyte Foodlyte, with its certified bio-compatibility. The standard EasyFerm Bio, with its HB glass, is directed at the Food & Beverage industry where CIP and SIP cycles occur frequently because the glass shows a very fast recovery leading to a shortened set-up time. The new versions with the PHI glass show an excellent chemical robustness at high pH values.

The LEVP (LE = Liquid Earth) versions have a stabilized sensor signal and an extended sensor diagnosis.



Benefits

- ➤ Specifically designed for sterile applications in Pharma and Biotechnology (Biocompatibility)
- ► Highly reliable measurements after steam sterilization, autoclavation and CIP cleanings
- ➤ Drift free measurements
- Ceramic diaphragm is an improved barrier of the electrode


Typical applications

- ▶ Bioreactors
- Downstream processes
- ▶ Brewhouse
- Gelatine manufacturing









Ordering Information

243632					
	Code	pH glas	SS		
	1	PHI			
	2	HB (reco	ommended	pH glass	type)
		Code	Electric	cal Conne	ctor
		1	VP 🚱		
		2	S8 🚱		
		3	Arc		
		4	Memose	ens 😉	
		5	K8 🚱		
		6	LEVP (o	nly for 120	and 225 mm length) 69
			Code	a-lengt	h (mm)
			1	120	
			2	160	
			3	200	
			4	225	
			5	325	
			7	425	
				Code	Temperature sensor
				1	Pt100 (VP, LEVP) (not applicable for Arc)
				2	Pt1000 (VP, LEVP) (not applicable for Arc)
	+	+	+	3	none (S8, K8) or given (Memosens, Arc)
243632 –					← Order Code



Accessories


pH buffers see page **▶** 100

Cables see page ▶ 108

Housings see page ▶ 126

MecoTrode





For more specifications see www.hamiltoncompany.com

The maintenance free MecoTrode sensor is designed for processes in the chemical industry with extreme pH values. The H glass type membrane glass provides a low alkaline error and stable measurement even at high temperatures.

Three high-performance ceramic diaphragms reduce the effect of flow potential in pipe mounting.



Benefits

- ➤ 3 high performance ceramic diaphragms for reduced flow potentials when mounted in pipes
- «H» glass for most accurate readings at high pH values or high temperatures
- ➤ Very good precision at low pH values (pH < 2)

Typical applications

- Water and Wastewater
- Industrial processes

Ordering Information a-length S8 VP 6 MS Arc MecoTrode H 120 238801 238437 242837 10110152* MecoTrode HF 120 242839

*Not for explosive environments

225

Accessories



pH buffers see page ▶ 100

Cables see page ▶ 108

Housings see page ▶ 126

242840





SENSORS **SENSORS**

OneFerm pH new



OneFerm pH VP 70

The OneFerm family of pH sensors is designed for applications in the single-use (SU) Pharmaceutical and Biotechnology Industries. Hamilton OneFerm sensors are the next step in the evolution of singleuse measurement. Their design solves some of the issues that commonly occur with reusable pH sensors that are inserted into the bag.

Specifically, Hamilton's single-use sensors combine the reliability and measurement stability of our longterm proven conventional sensors with the ease of use as an integral part of the bioreactor. The sensors retain the high accuracy performance even after gamma irradiation and a sufficient shelf life making it the ideal single-use solution.



Did you know... that with the reusable Arc Module SU pH a very stable digital signal can be achieved?

Benefits

- ► Specially designed for sterile application in SU Pharma
- ► Highly reliable measurements after gamma sterilization and dry storage even after short wet-in time (<30 min)
- ► Very low drift (<0.1 pH per week)
- ► Biocompatible materials (ISO 10993-5 and USP <87>)

Typical applications



Ordering Information





Specifications

Measuring range Process temperature

Pressure range (relative to ambient) Hygienic aspects

Diaphragm

O-ring



3 to 10 pH

4 to 50 °C 0 to 1 bar

pH-port)

Silicone

For more specifications see www.hamiltoncompany.com

HP Coatramic

Gamma irradiation up to 45 kGy (for the OneFerm sensors and the



	a-length	VP 6 / Pt100	VP 6 / Pt1000	VP 6 / NTC22	K8
OneFerm pH*	70	243216	243266	243235	_
	120	243217	243267	243236	243271
	160	10064894	10108674	10065001	10106075
	225	243218	243268	243237	243272
	325	243219	243269	243238	243273
	425	10101065	10089592	243239	243274

*Only for OEM integration available







Accessories



Cables see page ▶ 108

Of

SENSORS **SENSORS**

ChemoTrode / P ChemoTrode Bridge



The ChemoTrode is the most robust sensor to measure pH in demanding applications in pharmaceutical and chemical industries.


The ChemoTrode has a refill hole which allows refilling of the electrolyte and pressurization of the reference system. Its Everef-F reference cartridge ensures that the reference electrolyte remains free of silver and precipitation of proteins.

Did you know... that the ChemoTrode Bridge has an extended life time due to its special reference system?

Benefits

- Liquid electrolyte ensures fast response time and high
- ► Longer lifetime thanks to refillable electrolyte
- ► Everef-F reference cartridge extends electrode life in aggressive media

Typical applications



For more specifications see www.hamiltoncompany.com







	a-length	S7	VP 6 / Pt1000	VP 6 / Pt100	
ChemoTrode	120	238760	242700	_	
	150	238762	242701	_	
	200	238764	-	_	
	250	238766	242703	10069903	
ChemoTrode P	120	238761	243252	-	
	150	238763	243253	_	
	250	238767	243254	-	
ChemoTrode Bridge	120	238770	-	-	
(Non Ex)	150	238772	-	_	
	250	238776	_	_	

Accessories



pH buffers see page **≥** 100 Cables see page ▶ 108 **Housings** see page **▶** 126

- ► Industrial processes



29 28

Ordering Information



FermoTrode



) to 130 °C) to 4 bar
) to 4 bar
SIP
PHI
Skylyte
Everef-F

For more specifications see www.hamiltoncompany.com

The maintenance free FermoTrode sensors are designed for measuring pH in pharmaceutical and biotechnological industries and fit in the MasterFit and RetractoMaster housings. The Everef-F reference cartridge ensures that the reference electrolyte Skylyte remains free of silver and precipitation, and withstands steam sterilization.

It is not suited for contact with caustic soda like in CIP-cleanings or for use in media containing citric acid.

Benefits

- ➤ No air pressure required, no risk of empty reference electrolyte compartment
- ▶ 3 Coatramic diaphragms prevent clogging due to proteins
- Very long lifetime, stable calibration after sterilization and practically drift-free signals

Typical applications

- Biotechnology
- Pharmaceutical Industry



	a-length	S7	
FermoTrode	120	238480	
	150	238482	
	200	238484	
	250	238486	

Accessories



pH buffers see page ▶ 100

Cables see page ▶ 108

Housings see page ▶ 126

IEC TECEX

30 31

Ordering Information

IonoTrode



Specifications Measuring range 0 to 14 pH Process temperature -10 to 40 °C 0 to 0.5 bar or higher if Pressure range (relative to ambient) pressurization by side-arm pH glass 3 M KCI Electrolyte Reference system Everef Diaphragm Sleeve **EPDM** O-ring

For more specifications see www.hamiltoncompany.com

The lonoTrode sensor is designed for applications in ion weak media. The F glass membrane has a very low resistance, therefore the sensor can be used in samples with low conductivity, where it offers highest accuracy over a long period of time.

If there is a storage container with 3 M KCl attached via a tube to the side-arm of the lonoTrode, the flow-out of the electrolyte can be controlled with the sleeve diaphragm.



Did you know...

that the IonoTrode is designed for ion weak media with a low conductivity of only 0.2 µS/cm?

Benefits

- ► Offers highest accuracy over a long period of time
- ► Stable measurements in samples with low conductivity of at least 0.2 µS/cm
- ► Removable PTFE sleeve diaphragm to check electrolyte
- ➤ Side-arm attachment via tube to storage vessel containing 3 M KCl, and control of electrolyte flow with PTFE diaphragm ring

Typical applications

- Drinking Water Plants
- ▶ Boiler Feed Water

Ordering Information



	a-length	S 7
lonoTrode	120	238525

Accessories



pH buffers see page ▶ 100

Cables see page ▶ 108

Housings see page ▶ 126

SENSORS SENSORS

InchTrode



Specifications Measuring range 0 to 14 pH Process temperature -10 to 130 °C (flat membrane) 0 to 130 °C (cylindrical membrane) 0 to 10 bar (25 °C) Pressure range (relative to ambient) 0 to 6 bar (130 °C) pH glass HF (flat membrane) PHI (cylindrical membrane) Electrolyte Polisolve Reference system Everef-L Diaphragm Single Pore Temperature sensor Pt1000 in VP version Pt100 in fix cable version

For more specifications see www.hamiltoncompany.com

The InchTrode sensors are designed to measure pH in demanding applications in the paper making as well as in the chemical industries. The Single Pore liquid junction guarantees the best and fast measuring results because of direct contact between the sample and the Polisolve electrolyte.

The InchTrode sensors are easy to install without additional housing and have a robust PEEK shaft.

Did you know... that the InchTrode is available in two different sizes and with different membrane shapes?

Benefits

- ➤ Single Pore for direct sample contact with Polisolve electrolyte no clogging
- ➤ Very long-lasting reference system
- ► Robust PEEK shaft
- ► Simple installation without additional housing

Typical applications

- Pulp and Paper industry
- ▶ Water and Wastewater

Ordering Information



F = Flat membrane P = Cylindrical membrane

 $C = Fix \ cable$

Accessories



pH buffers see page ▶ 100


Cables see page ▶ 108

Housings see page ▶ 126

IEC IECEX



Polilyte Pro Polyplast Pro



Specifications Measuring range 0 to 14 pH Process temperature Polilyte Pro: -10 to 60 °C Polyplast Pro: -10 to 40 °C 0 to 6 bar Pressure range (relative to ambient) pH glass Polilyte Pro: HF Polyplast Pro: V Polisolve Electrolyte Reference system Polilyte Pro: Everef-B Polyplast Pro: Ag/AgCl Diaphragm Single Pore Pt1000 in VP version Temperature sensor O-ring Polilyte Pro: EPDM Polyplast Pro: EPDM

For more specifications see www.hamiltoncompany.com

The maintenance free Polilyte Pro and Polyplast Pro sensors are designed for pH measurement in water applications, especially in low conductivity samples, e.g. wastewater, fish farming, ground water, etc.

The Single Pore liquid junction guarantees best measurement results because of direct contact between the sample and the Polisolve electrolyte – clogging is nearly impossible. The Polyplast Pro sensor comes with a robust plastic shaft and glass bulb protection.

Did you know...

that the Polilyte Pro has the
HF resistant pH glass?

Benefits

- ➤ Single Pore for direct sample contact with Polisolve electrolyte
- ▶ No clogging
- ► Fast response even in low conductivity media
- ► Easy maintenance due to non-refillable electrolyte

Typical applications

- Wastewater applications
- Fish farming
- Ground water



	a-length	S 8	VP 6
Polilyte Pro	120	238411	238417
Polyplast Pro	120	238408	-

Accessories



pH buffers see page ▶ 100

Cables see page ▶ 108

Housings see page ▶ 126

IEC IECEX



36 37

Ordering Information

SENSORS SENSORS

Liq-Glass PG EasyControl



The maintenance free Liq-Glass PG and the EasyControl sensors are entry level sensors for chemical or waste water applications and low process temperatures. They show good behaviour in samples with low conductivity.

Did you know...

that the Easy Control is also
available as ORP sensor?

Benefits

- ➤ Suitable for low conductivity media
- ► Easy maintenance due to non-refillable electrolyte
- ► Liq-Glass PG has 3 ceramic diaphragms for reduced flow potentials

Typical applications

- Wastewater applications
- ▶ Fish farming
- Ground water
- Swimming Pools

Ordering Information



Accessories



pH buffers see page ▶ 100

Cables see page ▶ 108

Housings see page ▶ 126

IEC TEĈE



Specifications Liq-Glass PG: 1 to 12 pH Measuring range EasyControl: 0 to 14 pH Lig-Glass PG: -5 to 60 °C Process temperature EasyControl: 0 to 60 °C 0 to 2 bar Pressure range (relative to ambient) Liq-Glass PG: F pH glass EasyControl: HF Electrolyte Liq-Glass PG: Viscous 3 M KCI-LR EasyControl: Gel electrolyte Lig-Glass PG: Everef Reference system EasyControl: Ag/AgCl Diaphragm Ceramic Liq-Glass: EPDM O-ring EasyControl: EPDM

For more specifications see www.hamiltoncompany.com



ORP (Oxidation Reduction Potential) is a common measurement in biochemistry, environmental chemistry and water quality. In the biochemical perspective, an oxidizing chemical pulls electrons away from the cell membrane which means it can be destabilized and leaky. The rapid death of a cell is the consequence of a destroyed membrane. The ORPs of natural systems like aerated surface water, rivers, lakes, rainwater and acid mine water usually have oxidizing conditions leading to positive potentials. Submerged soils, swamps and marine sediments, where air supply has its limitations, reducing conditions are the norm leading to negative potentials. For water system monitoring, the ORP value provides the operator with a rapid and single-value assessment of the disinfection potential of water in the postharvest system. This enables the operator to assess the activity of the applied disinfectant rather than the applied dose.


ORPs in aqueous solutions are determined by measuring the potential difference between an inert sensing electrode in contact with the solution and a stable reference electrode. The reference electrode is connected to the solution by a salt bridge. It has a known potential and is made of silver chloride or saturate calomel. Platinum is frequently used for the sensing electrode.

The Oxygen-Reduction Potential, also known as Redox Potential describes the tendency of a chemical species or a solution to acquire electrons and therefore to be reduced. Each species has its own reduction potential. It is measured in Volts (V) or mV.



Polilyte Plus ORP





SENSORS

For more specifications see www.hamiltoncompany.com

The maintenance free Polilyte Plus ORP sensors are designed to withstand demanding applications in chemical and petrochemical industries. Monitoring the ORP value is becoming increasingly important in many applications, especially harsh chemical environments or high alkaline wastewater. Because of its Single Pore diaphragms you will never have liquid junction problems and total breakdowns. The Polilyte Plus ORP sensors demonstrate reliable reproducible measurement accuracy in highly alkaline solutions as well as in samples with low conductivity. Additionally, the Everef-L reference cartridge ensures a long lifetime.

Benefits

- ▶ 2 Single Pores prevent clogging and ensure reliable measurements
- Minimal diffusion potenital
- ► Highly reproducible measurements and very stable over a long period of time
- ► Resistant against solvents, strong acids and bases

Typical applications







ORP buffers see page **▶** 101

Cables see page ≥ 108

Arc Accessories see page **≥** 117

Housings see page **▶** 126







EasyFerm Plus ORP



Specifications Measuring range ± 2000 mV (Arc: ± 1500 mV) Process temperature 0 to 140 °C (Arc: analog 0 to 110 °C, digital 0 to 140 °C) Pressure range 0 to 6 bar (relative to ambient) Hygienic aspects Autoclavable, CIP, SIP ORP element Pt ring Electrolyte Phermlyte Everef-F Reference system Diaphragm **HP** Coatramic O-ring **EPDM**

For more specifications see www.hamiltoncompany.com

The EasyFerm Plus ORP sensors are designed to withstand demanding applications in pharmaceutical and chemical industries. It is supplied with a prepressurized electrolyte which prevents the diffusion of sample into the sensors. The Everef-F reference cartridge ensures that the Phermlyte reference electrolyte remains free of silver and precipitation.

Measuring the ORP value is getting more and more important in the branches mentioned above.

Benefits

- ► Pre-pressurized reference electrolyte ensures a clog-free diaphragm
- ► Almost drift-free measurement
- ► Stable measurement signals after steam sterilization, autoclavation and CIP cleanings
- ► Large platinum ring

Typical applications

- Bioreactors
- ► Industrial processes
- Downstream processes

Ordering Information





	a-length	\$8	Arc
EasyFerm Plus ORP	120	243187	243050
	225	243188	243051
	325	-	243052
	425	_	243053

Accessories



ORP buffers see page **▶** 101

Cables see page ▶ 108

Arc Accessories see page ▶ 117

Housings see page **▶** 126









ChemoTrode ORP



Measuring range	± 2000 mV
Process temperature	0 to 130 °C
Pressure range (relative to ambient)	0 to 6 bar
ORP element	Pt ring
Electrolyte	Viscous 3 M KCI-LR
Reference system	Everef-F
Diaphragm	HP Ceramic

For more specifications see www.hamiltoncompany.com

The ChemoTrode ORP is the most robust sensor to measure the oxidation-reduction potential in demanding applications in pharmaceutical and chemical industries. The ChemoTrode ORP has a refill hole which allows refilling the electrolyte and pressurization of the reference electrolyte. Its Everef-F reference cartridge ensures that the reference electrolyte remains free of silver and precipitation of proteins.

Benefits

- ► Liquid electrolyte ensures fast response time and
- ► Longer lifetime thanks to refillable electrolyte
- ► Everef-F reference cartridge extends electrode life in aggressive media

Typical applications

- Industrial processes
- Mining Industry
- Pulp and Paper industry
- ► Fermentations

Ordering Information



Accessories



ORP buffers see page **▶** 101 **Cables** see page **▶** 108

Housings see page **≥** 126







OxyTrode Pt



Specifications Measuring range ± 2000 mV 0 to 130 °C Process temperature 0 to 16 bar (25 °C) Pressure range 0 to 6 bar (130 °C) (relative to ambient) ORP element Pt wire Viscous 3 M KCI-Pharma, blue Electrolyte Reference system Everef HP ceramic Diaphragm **EPDM** O-ring

For more specifications see www.hamiltoncompany.com

The maintenance free OxyTrode Pt is an ORP sensor designed for processes in the chemical industry and for applications in wastewater treatment. Three high-performance ceramic diaphragms reduce the effect of flow potential in pipe mounting.



Benefits

- ➤ 3 high performance ceramic diaphragms for reduced flow potentials when mounted in pipes
- ▶ Platinum wire coil welded onto the glass

Typical applications

- Water and Wastewater
- Industrial processes

Ordering Information



Accessories



ORP buffers see page ≥ 101
Cables see page ≥ 108

Housings see page **▶** 126







Polilyte RX Polyplast Pro RX



Specifications ± 2000 mV Measuring range Process temperature Polilyte Pro: -10 to 60 °C Polyplast Pro: -10 to 40 °C 0 to 6 bar Pressure range (relative to ambient) ORP element Pt-wire Polisolve Electrolyte Reference system Polilyte Pro: Everef-B Polyplast Pro: Ag/AgCl Diaphragm Single Pore O-ring Polilyte RX: EPDM Polyplast Pro RX: EPDM

For more specifications see www.hamiltoncompany.com

The maintenance free Polilyte RX and Polyplast Pro RX sensors are designed for ORP measurement in water applications and low conductivity samples, e.g. wastewater, fish farming, ground water, etc.

The Single Pore liquid junction guarantees best measurement results because of direct contact between the sample and the Polisolve electrolyte – clogging is nearly impossible. The Polyplast Pro sensor comes with a robust plastic shaft and glass bulb protection, making it one of our most economical and longest lasting sensors.

Benefits

- ➤ Single Pore for direct sample contact with Polisolve electrolyte
- ▶ No clogging
- ► Fast response even in low conductivity media
- ► Easy maintenance due to non refillable electrolyte

Typical applications

- Wastewater applications
- Fish farming
- ► Ground water

Ordering Information



Accessories



ORP buffers see page ▶ 101
Cables see page ▶ 108

Housings see page **▶** 126







EasyControl ORP



Specifications Measuring range \pm 2000 mV Process temperature 0 to 60 °C 0 to 2 bar Pressure range (relative to ambient) ORP element Pt-wire Electrolyte Gel electrolyte Ag/AgCI Reference system Diaphragm Ceramic **EPDM** O-ring

For more specifications see www.hamiltoncompany.com

The maintenance free EasyControl ORP is an entry level ORP sensor for chemical or wastewater applications and low process temperatures.

It is also often used in swimming pools to control the disinfection with chlorine. They show also good behavior in samples containing few ions, with respectively low conductivity.

Benefits

- ➤ Suitable for low conductivity media
- ► Easy maintenance due to non refillable electrolyte

Typical applications

- Wastewater applications
- ▶ Fish farming
- Ground water
- Swimming Pools

Ordering Information



Accessories



ORP buffers see page **▶** 101 **Cables** see page **▶** 108


Housings see page **≥** 126



The electrical conductivity is important for the characterization of liquids in different kinds of processes. In aqueous solutions the conductivity is caused by the decomposition of dissolved acids, bases or salts into positive cations and negative anions. In ultra-pure water, where no ions, except very few $\rm H_3O^+$ and $\rm OH^-$, are present, the conductivity is extremely low. This intrinsic conductivity of water represents the lower border of the conductivity scale.

The electrical conductivity is determined by a resistivity measurement when an alternating voltage is applied to a measurement cell that consists of two or four electrodes. To compensate for the geometry of the conductivity cell a cell constant is used. This constant is either known or determined by means of conductivity standards.

Electrical conductivity is the reciprocal of electrical resistivity, and measures a material's ability to conduct an electric current. Its SI unit is Siemens per meter (S/m). For the measurement of the conductivity of a solution it's common to use μ S/cm or mS/cm.



Cond SENSORS

SENSORS

Conduce 4UxF [family]



Specifications 1 μ S/cm to 300 mS/cm Measuring range Measurement Principle 4 pole contacting Process temperature -20 to 150 °C (Arc: analog 0 to 110 °C, digital 0 to 140 °C) 0 to 20 bar (135 °C) Pressure range (relative to ambient) 0 to 10 bar (150 °C) Autoclavable, CIP, SIP Hygienic aspects Cell constant 0.36/cm Material of electrodes S = Stainless steel 1.4435 H = Hastelloy C 2.4602 T = TitaniumPt = Platinum O-ring EPDM (other versions available on request)

For more specifications see www.hamiltoncompany.com

The Conducell 4UxF sensors are suited for measurements in hygienic applications. All wetted parts are FDA-approved, can be cleaned easily and withstand CIP cleanings and autoclavations. The sensors show a very good linearity over a broad measuring range.

They are available with different process connections such as Varivent®.

The Conducell 4USF with stainless steel electrodes is most common. This sensor is suitable for various applications in biopharma, water or food industry. The newly implemented lengths are perfectly designed for flow cells e.g. in downstream applications.

All plastic materials are compliant with the order FU 10/2011.

Benefits

- ► Very good linearity, especially for applications with sharp variations in conductivity
- ► All wetted parts are FDA-compliant
- ► Sensor is very easy to clean due to the forward facing, flush arrangement of electrodes
- ► Specifically designed for sterile applications in Pharma and Biotechnology

Typical applications

Ordering Information

243590					
***************************************	Code	Electro	de Materi	al	
	1	Stainless	Steel 1.4	435	
	2	Platinum	(not for Tr	riclamp)	
	3	Stainless	Steel 2.4	602	
	4	Titanium	(not for Tr	iclamp)	
		Code	Electric	al Conne	ector
		1	Arc		
		2	VP 😉		
			Code	a-lengt	th (mm)
			1	120 (PG	
			2	225 (PG	313,5)
			3	325 (PG	313,5)
			4	425 (PG	h13,5)
			5	30 (PG1	(3,5)
			6	60 (PG1	(3,5)
		7		21 - Triclamp 1.5"	
				Code	O-ring Material
	+	+	+	1	EPDM
243590 -					← Order Code



Accessories



Conductivity Standards see page **≥** 102 Cables see page ≥ 108 **Housings** see page **▶** 126





sensors sensors

Conducell SU new



Hamilton's single-use conductivity monitoring system is comprised of the reusable Arc Module Cond-P SU and a single-use sensor patch Conducell-P SU. The Conducell-P SU is integrated within the single-use container by the container manufacturer.

Unlike other single-use conductivity solutions, Hamilton's reusable Arc Module enables a compact and cost-effective measurement solution without sacrificing accuracy or precision. A standard measuring loop consists of a sensor element (Conducell-P SU), which is connected directly to the electronic (Arc Module Cond-P SU) to enable disturbance free measurement signals.



Did you know... that with the reuseable Arc Module and the precalibrated sensor a ready to use system can be achieved?

Benefits

- ➤ Specially designed for sterile application in SU Pharma and Biotechnology
- ► Highly reliable measurements after gamma sterilization and dry storage even after short wet-in time (<30 min)
- ► Biocompatible materials

Typical applications

 Mixing bags for buffer preparation, virus inactivation or intermediate storage





0.1 to 300 mS/cm

(for the disposables)

Gamma irradiation up to 50 kGy

4 to 50 °C

0 to 1 bar

1.31/cm

Pt = Platinum

Arc Module Cond-P SU Conducell-P SU*
10071707 10076677

Specifications

Measuring range

Process temperature
Pressure range

(relative to ambient)

Hygienic aspects

Material of electrodes

For more specifications see www.hamiltoncompany.com

Cell constant

Measurement Principle 4 pole contacting

*Only for OEM integration available



Accessories

Conductivity Standards see page ▶ 102 Cables see page ▶ 108

Of

Cond

Conducell 4US



The Conducell 4US 4-pole conductivity sensors are designed for different process connections such as Triclamp or G 11/4" with various O-ring positions.

The sensors show a very good linearity over a broad range of conductivities.

The Conducell 4US 4-pole sensor can easily bei cleaned and is suitable for steam sterilization, autoclavation and CIP cleanings.


All plastic materials are compliant with the order EU 10/2011.

Benefits

- ► Very good linearity, especially for applications with wide variations in conductivity
- ► All wetted parts are FDA-compliant
- ➤ Sanitary: Sensor is easy to clean
- ► O-ring position can be chosen individually

Typical applications

- Fermentation
- Chemical industry



Specifications	
Measuring range	0.1 µS/cm to 500 mS/cm
Measurement Principle	4 pole contacting
O-ring position	22 to 55 mm
Process temperature	-20 to 135 °C
Pressure range (relative to ambient)	0 to 6 bar
Hygienic aspects	CIP, SIP
Cell constant	0.147/cm
Material of electrodes	Stainless steel 1.4435
O-ring	EPDM (other versions available on request)

For more specifications see www.hamiltoncompany.com

Ordering Information



	a-length	5 m fix cable
Conducell 4US-G125	variable	237700-OP
Conducell 4US-T150-50	50	237750
Conducell 4US-T150-100	100	237760

Accessories



Flow-through cell PEEK TC 1.5" Ref 237931
 This flow through cell made of FDA approved PEEK facilitates insertion of Conducell 4US-T150-50 in pipework.

Conductivity Standards see page ▶ 102 Safety Socket see page ▶ 152





sensors

Conducell UPW



Specifications Measuring range 0.01 to 1500 µS/cm Measurement Principle 2 pole contacting Arc: analog 0 to 110 °C, Process temperature digital 0 to 130 °C Pressure range 0 to 10 bar (130 °C) (relative to ambient) Autoclavable, CIP, SIP Hygienic aspects Cell constant < 0.1/cmMaterial of electrodes Stainless Steel DIN 1.4435 Surface quality $R_a < 0.4 \ \mu m \ (N5)$ O-ring EPDM (other versions available on request)

For more specifications see www.hamiltoncompany.com

The Conducell UPW 2-pole conductivity sensors are designed for the use in liquids with very low conductivity, i.e. Ultra Pure Water, Pure Water and Water for Injection, particularly in the pharmaceutical and chemical industry.

Conducell UPW sensors are available with different process connections such as TriClamp 1.5", PG 13.5.

All plastic materials are compliant with the order EU 10/2011.




Benefits

- ► Sanitary design: all wetted parts are FDA approved
- ► Easy cleanable
- ► Intelligence in the sensor: fully compensated measurement signals
- ► Easy handling due to user-friendly interface

Typical applications

- ➤ Ultra Pure Water
- Pure Wate
- Water for Injection

Ordering Information





Accessories



UPW Simulator Ref 243580

Traceable resistor to verify the Arc module acc. to USP <645>

Conductivity Standards see page **▶** 102

Cables see page ▶ 108

Arc Accessories see page **▶** 117

Housings see page **▶** 126



SENSORS SENSORS

Conducell 2DC-PG



Measuring range	10 μS/cm to 20 mS/cm
Measurement Principle	2 pole contacting
Process temperature	-5 to 80 °C
Pressure range (relative to ambient)	0 to 6 bar
Cell constant	1/cm
Material of electrodes	Graphite
O-ring	EPDM (other versions available on request)

For more specifications see www.hamiltoncompany.com

The Conducell 2DC sensor is constructed in a simple way and is best suited for measurements in clean solutions and non-critical applications. Contaminants, such as lime, will affect the measurement.





	a-length	5 m fix cable
Conducell 2DC-PG 120	120	237610

Benefits

- ▶ 2 large graphite electrodes for stable measurements
- ► Mechanically-stable plastic shaft
- ► Easily cleanable

Typical applications

Water and Wastewater

Accessories



Conductivity Standards see page ▶ 102
Housings see page ▶ 126



Biological processes are increasingly important in biotechnical and pharmaceutical industries. The variability of living organisms is often very high, making the culture process difficult to standardize. Extensive process optimization and control are required for stable cell cultures, fermentations and improved yield. Today bioprocess development relies on labor intensive sampling and offline measurements that also lack the necessary granularity to fully optimize the yield. The available on-line measurements of pH and dissolved oxygen are not linked to the cell status

On-line monitoring of cell density provides the continuous information necessary to optimize control and yield beyond what is possible off-line. Hamilton now offers sensors for continuous cell density measurement. The Incyte Arc permittivity sensor delivers information on viable cell density whereas the Dencytee sensor measures total cell density via turbidity. In combination with our advanced Arc pH and dissolved oxygen probes, permittivity and turbidity sensors provide all relevant information on the process of mammalian, yeast and high density bacteria cultures. This enables better

SENSORS **SENSORS**



Incyte Arc



Viable Cell Density Measurements in real-time provide actionable data to automate process control. This is not possible by relying on off-line measurements, which only provide a window into the past.

Incyte Arc is Hamilton's next-generation viable cell density sensor, offering high-fidelity permittivity measurements comes now paired with integrated microtransmitters that leverage ArcAir technology. Arc Wi 2G Adapter BT (REF 243470) is required to output an analog 4-20 mA signal from the digital Modbus communication. Arc Wi 1G Adapter BT (REF 242360) is required with Arc Modbus OPC Converter (REF 10089359) to enable an OPC communication.



Did you know... Incyte Arc is now part of the Hamilton Arc family providing a digital Arc Modbus signal directly

Benefits


- Specific for viable cells
- ➤ Suitable for cell culture
- ► Insensitive to microcarriers and cell debris
- ► No more hidden events
- ► Optimization of feeding strategy and yield

Typical applications



from the sensor? 99

Ordering Information



Measuring Range 5 x 10⁵ to 8 x 10⁹ cells/mL (Mammalian) Conductivity range 0.5 to 80 mS/cm Measuring principle Permittivity Process temperature 0 to 60 °C Pressure range 0 to 12 bar Hygienic aspects Autoclavable, CIP, SIP

EPDM

For more specifications see www.hamiltoncompany.com

Specifications

O-ring





Conductivity standard 12880 µS/cm, Basic Line Ref 238988

Solution B for Incyte Verification

Ref 243742

Cables see page ▶ 108

Arc Accessories see page **≥** 117

Housings see page **≥** 126



SENSORS **SENSORS**



Specifications Measuring Range

Conductivity range Measuring principle

Process temperature

(relative to ambient) Hygienic aspects

Material of electrodes

For more specifications see www.hamiltoncompany.com

Pressure range

Incyte SU new

Accurately analyzing the characteristics of viable cells during bioprocess is crucial. Only viable and healthy cells are producing the product of interest. Today these characteristics are monitored by labor intensive offline samplings.

Analyzing cell characteristics online provides deep insight into the bioprocess. It allows stable process control, fast optimization and reduces the risk of sampling errors. The Incyte SU sensor is especially designed for measuring viable cells during mammalian cell culture, yeast and high density bacterial fermentation.

The measurement principle of Incyte sensors is based on permittivity. Viable cells behave like little capacitors and their polarization and depolarization in an alternating electrical field is measured. This signal can be correlated to the viable cell density. This method is insensitive to cell debris and microcarriers because only viable cells can be polarized.

A measuring Unit consists of an sensor element (Incyte-P SU) and an electronic (Arc Module Incyte-P SU), which converts the analog measurement to a stable digital signal.



Did you know... that Hamilton is the only provider of all relevant parameters in single use and re-usable technology for cell culture & fermentations: viable cell density, pH and DO?

Benefits

- ➤ Specific for viable cells
- ► Suitable for cell culture and fermentation
- ► Insensitive to microcarriers and cell debris
- ➤ No more hidden events
- ► Optimization of feeding strategy and yield
- ► Biocompatible materials

Typical applications





*Only for OEM integration available

5 x 10⁵ to 8 x 10⁹ cells/mL

Gamma irradiation up to 50 kGy

(for the disposables)

(Mammalian) 1 to 50 mS/cm

Permittivity

4 to 50 °C

0 to 1 bar

Platinum

Accessories



Conductivity Standards see page **▶** 102 Cables see page ≥ 108



SENSORS **SENSORS**



Dencytee Arc new



Specifications e.g. 0 to 200g/I cell dry weight yeast Measuring Range 0 to 4 AU 0 to 30'000 NTU Measuring principle Transmission and Reflection (incl. temperature compensation, daylight filter and subtraction) Wavelength 860 nm

Dencytee Arc sensors perform on-line measurement of Total Cell Density in solution. The sensor is designed to measure the turbidity of the cell suspension. The measurement is made at NIR (near-infrared) wavelengths so it is insensitive to changes in media color.

All particles and molecules that scatter light at 860 nm will be detected, including living and dead cells as well as cell debris. The sensor is also very effective after inoculation when cells are expanding quickly but concentrations are low, making capacitance-based readings less reliable.

Dencytee Arc sensors provide a robust connection directly to the Process Control System without the need for an additional external transmitter.



Did you know... to be able to measure low and high cell density at a high quality signal the sensor is able to measure the transmitted as well as the reflected light of the cells. 99

Benefits

- ► Simple On-Line Measurement of Cell Growth
- ► Reliable Values During the Growth Phase
- ► Improved Linearity at High Concentrations
- ► Early Detection of Process Deviations

Typical applications



Ordering Information

	a-length	Arc	
Dencytee RS485	120	10064919-11	
Dencytee RS485	225	10064919-12	
Dencytee RS485	325	10064919-13	
Dencytee RS485	425	10064919-14	

Accessories



Dencytee Maintenance Tool Kit Ref 10146924

Cables see page ▶ 108

Arc Accessories see page **▶** 117

Housings see page **▶** 126



75



Dissolved carbon dioxide (DCO₂) is a critical process parameter (CPP) in biopharma production processes according to PAT guidelines. By influencing other parameters such as extracellular and intracellular pH, it has an effect on different metabolic pathways which are involved in cell growth or in product formation and quality.

In the past, continuous in-line monitoring of DCO₂ has only been possible through electrochemical sensors that are based on the Severinghaus principle and measure the DCO₂ concentration indirectly. The result is significant maintenance effort and multiple sources of drift that must be compensated by time-consuming product calibration.

Now, Hamilton has introduced a completely new way to measure DCO₂: The new in-line sensor CO₂NTROL is a maintenance free, solid-state sensor that directly measures DCO₂ resulting in better measurement accuracy and lower cost of ownership.

sensors sensors

CO₂NTROL new



The Solid State Sensor directly measures DCO₂ and provides maintenance free, real-time, and in-line control of this new critical process parameter.

Unlike traditional sensors that are based on the electrochemical Severinghaus principle, CO_2NTROL is a pure direct measurement in a solid state design: CO_2 molecules diffuse into a gas permeable membrane where the sensor measures the absorption of CO_2 -specific Mid-IR wavelengths. This absorption correlates to the partial pressure of CO_2 in the media.

CO₂NTROL's hygienic design makes it compliant with requirements of biopharma applications. The sensor is EHEDG approved (EL Class I, test executed with Hamilton hygienic socket REF 242545) and is ready for GMP compliance. Embedded electronics convert the MIR CO₂ measurement into standard digital and analog signals that are easily integrated into your control strategy.

Arc Wi 2G Adapter BT (REF 243470) is required to output an analog 4-20 mA signal from the digital Modbus communication.



Did you know...

Hamilton is the first and only supplier to bring the maintenance-free optical IR technology into a SIP/CIP compliant 12mm CO₂ sensor 99

Benefits

- ► Maintenance-free
- ► Simple calibration
- ► Hygienic design: SIP/CIP compatible, autoclavable
- Inverted installation possible
- ▶ Direct measurement of CO₂ no ammonia interference





Typical applications

► Biopharma Cell Cultures and Fermentations





*CO₂NTROL 225 have, in reality, a shaft length of 215 mm. This ensures optimal rinsing in replaceable armatures, such as Retractex.

Accessories



Calibration Station Ref 243575

Cables see page ▶ 108

Arc Accessories see page **▶** 117


Housings see page **▶** 126

Specifications Measurement Principle Optical – CO₂ Absorption in Middle Infrared (MIR) 5 to 1000 mbar Measuring Range or 0.5 to 100 %-Vol or 7.5 to 1500 mg/L (in liquid phase at 101.3 kPa and 25 °C) Diameter 12 mm **Process Connection** PG 13.5 Wetted Parts Stainless Steel 1.4435, EPDM (Ethylene propylene elastomer), FDA compliant silicone Surface Quality $R_a < 0.4 \mu m (N5)$ Steam Sterilizable Yes Autoclavable Yes CIP Yes Operating -10 to 60 °C temperature range



The partial pressure of dissolved oxygen (DO) plays an important role in many biological, chemical and physical processes. The amount of dissolved oxygen is also important for the safety and the quality of many other industrial processes.

The most common technologies to measure DO are the classical amperometric and the modern optical method. Classical amperometric Clark cells, where cathode and anode are separated from the sample by a gas permeable membrane, generate an electrical current proportional to the oxygen partial pressure of dissolved oxygen. The oxygen is reduced in the sensor, catalyzed by an electrolyte at a platinum cathode. At the anode silver is oxidized. In contrast to the Clark cells the optical measurement is based on the luminescence of a luminophore that absorbs photons and releases a part of the absorbed energy by emission of photons with a higher wavelength. Oxygen quenches this process by transferring the energy partially by collision. The more oxygen present the more quenching is observed. Hamilton measures the phase shift between excitation and emission across a population of light pulses in order to achieve the highest accuracy and widest operating range. The difference in the intensity of both waves is used for online sensor diagnostics.



SENSORS **SENSORS**

VisiFerm RS485 family new







The VisiFerm RS485 is the first optical oxygen sensor with integrated opto-electronics, having the full functionality of a measuring device with self-diagnostics. It is steam sterilizable, autoclavable and CIP compatible. The VisiFerm requires less maintenance than a classical oxygen sensor as it does not have a mechanically sensitive membrane or a corrosive electrolyte.

Did you know... that Hamilton invented the first optical DO sensor in 12 mm format?

Benefits

- ► Reliable and robust optical measurement
- ► No fragile membrane with a solid sensor cap
- ► No polarization time required
- Instantly stable values, low drift, quick response
- ► Electrolyte-free, so no leakage
- ► Convenient precalibration in the laboratory, because data is stored in the sensor head
- ► Calibration, verification, and maintenance data accessible via ArcAir app

Typical applications

- ▶ Biotechnical fermentation

Ordering Information



ODO Cap H0 + H3: For general application in biotechnology, water treatment and monitoring as well as in breweries, wineries and soft drink processing.

ODO Cap H2 + H4: Designed for fermentation processes where sterilization in place (SIP) is performed in media containing higher amounts of lipophilic compounds. It comes with a hygienic design.

Accessories



• ODO Cap H0 Kit Ref 243515

Specifications

Measuring range Measurement Principle

Response time t98% Process temperature

Operating voltage

(relative to ambient) Hygienic aspects

Pressure range

Surface Quality

Material

O-ring

4 ppb to 25 ppm (DO)

10 to 27 VDC max. 1.5W

Autoclavable, CIP, SIP

Stainless steel 1.4435

 $R_a < 0.4 \mu m (N5)$

-1 to 12 bar

EPDM

For more specifications see www.hamiltoncompany.com

< 30 s at 25 °C, from air to nitrogen

-20 to 140 °C, the sensor provides no DO reading above 85 °C

Oxygen dependent luminescence quenching

- **ODO Cap H2 Kit** Ref 243505
- ODO Cap H3 Kit Ref 10068400
- ODO Cap H4 Kit Ref 10078261

Cables see page ▶ 108

Housings see page ▶ 126

Arc Accessories see page ▶ 117

Of

SENSORS

VisiFerm DO SU new



Specifications Measuring range 4 ppb to 25 ppm (DO) Measurement Principle Oxygen dependent luminescence quenching Response time t98% < 30 s at 25 °C, from air to nitrogen Process temperature 4 to 50 °C Operating voltage 7 to 30 VDC max. 1 W Hygienic aspects Gamma irradiation up to 50 kGy (for the disposables) O-ring **EPDM**

For more specifications see www.hamiltoncompany.com

Hamilton's single-use dissolved oxygen monitoring system is comprised of the reusable VisiFerm DO SU and a single-use optical dissolved oxygen sensor cap. The cap is integrated with the single-use container by the container manufacturer.

Hamilton's reusable sensor element enables a compact and cost-effective measurement solution without sacrificing accuracy or precision. A standard measuring loop consists of a sensor element, which is connected to the VisiFerm DO SU.



Benefits

- ► Specially designed for sterile application in SU Pharma and Biotechnology
- ► Highly reliable measurements after gamma sterilization and dry storage even after short wet-in time (<30 min)
- ► Very low drift
- ► Biocompatible material

Typical applications

- ► SU bioreactors (bag application)
- SU bioreactors (rigid containers)
- SU mixer (fill and finish application)



*Only for OEM integration available

Accessories



Silicone Sleeve (for ODO Cap S3) Ref 10114324

.....

Cables see page ▶ 108

Arc Accessories see page **▶** 117



SENSORS

VisiFerm mA family



Specifications 4 ppb to 25 ppm (DO) Measuring range Measurement Principle Oxygen dependent luminescence quenching Response time t98% < 30 s at 25 °C, from air to nitrogen Process temperature -20 to 140 °C, the sensor provides no DO reading above 85 °C Operating voltage 18 to 30 VDC Pressure range -1 to 12 bar (relative to ambient) Hygienic aspects Autoclavable, CIP, SIP Surface Quality $R_a < 0.4 \mu m (N5)$ Stainless steel 1.4435 Material O-ring **EPDM**

For more specifications see www.hamiltoncompany.com

The VisiFerm mA is the optical dissolved oxygen (DO) sensor for use in explosive environment. VisiFerm mA optical technology improves the measuring performance and simplifies maintenance. Improvements compared to conventional electrochemical (amperometric) sensors include flow independence, rapid startup with no polarization time, and simplified maintenance.

Designed especially for production environments, the VisiFerm mA is a 2-wire sensor with 4-20 mA standard or digital HART signal output, and ATEX & IECEx approval. The VisiFerm mA mitigates the negative effects of aging, temperature, and photobleaching in order to reduce the frequency of calibration and deviation reports.

Did you know... that Hamilton invented the first optical DO sensor in 12 mm format?

Benefits

- ► Reliable and robust optical measurement in hazardous
- ► Longer cap and sensor life
- ► Less frequent calibrations
- ► Easy installation with 2-wire connection
- ▶ Direct analog 4-20 mA or digital HART communication
- ► Calibration, verification, and maintenance data accessible via ArcAir app

Typical applications

- Explosive atmospheres environment
- Fermentation
- Wort agration in broweries

Ordering Information

10070760					
•••••	Code	Interfac	ce		
	1	mA/HAF	RT		
		Code	a-lengt	h (mm)	
		1	120		
		2	160		
		3	225*		
		4	325		
		5	425		
			Code	ODO Ca	ар
			1	НЗ	
			2	H4	
				Code	Wetted Parts
	+	+	+	1	EPDM
10070760 -					← Order Code

*The VisiFerm mA 225 have, in reality, a shaft length of 215 mm. This ensures optimal rinsing in retractable armatures, such as Retracte

ODO Cap H3: For general application in biotechnology, water treatment and monitoring as well as in breweries, wineries and soft drink processing.

ODO Cap H4: The ODO Cap H4 is designed for fermentation processes where sterilization in place (SIP) is performed in media containing higher amounts of lipophilic compounds. It comes with a hygienic design.

Accessories



- **ODO Cap H3 Kit** Ref 10068400
- ODO Cap H4 Kit Ref 10078261

Cables see page ≥ 108

Housings see page ≥ 126

HART O







SENSORS **SENSORS**

VisiTrace RS485 mily new







VisiTrace RS485 120

The VisiTrace RS485 is designed to measure dissolved oxygen in the low ppb ranges in brewing applications, notably during filtration, and filling. In addition, the special designed ODO Cap L1 for breweries is stabilized against standard disinfectant solution with active chlorine and chlorine dioxide. This is powerful during measurements in breweries, which may not allow for calibration after every CIP.

With the transmitter integrated, the intelligent VisiTrace RS485 sensor provides more reliable measurements.



Did you know... that the VisiTrace RS485 is the only optical DO sensor that withstands chlorine and chlorine dioxide for a long time? 99

Benefits

- ► For measurements from 0 to 2000 ppb
- ► Stable against chlorine and chlorine dioxide
- ► Rapid start-up with no polarization
- ► Flow and CO₂ independent readings
- ► Robust design for high flow rates

Typical applications

Ordering Information

10140043									
	Code	Interfac	ce						
	1	RS485							
		Code	a-lengt	h (mm)					
		1	120						
		2	160						
		3	225*	225*					
		4	325						
		5	425						
			Code ODO Cap						
			1	L1					
				Code	Wetted Parts				
	+	+	+	1	EPDM				
10140043 -					← Order Code				

Specifications

Measuring range Measurement Principle

Response time t_{98%}

Operating voltage

(relative to ambient) Hygienic aspects

Pressure range

Surface Quality

Material

O-ring

Process temperature

0 to 2000 ppb (DO)

quenching

-1 to 12 bar

EPDM

For more specifications see www.hamiltoncompany.com

 $R_a < 0.4 \mu m (N5)$

Autoclavable, CIP, SIP

Stainless steel 1.4435

Oxygen dependent luminescence

-20 to 140 °C, the sensor provides no DO reading above 85 °C 10 to 27 VDC max. 1.5W

< 20 s in gas; < 90 s in water

*The VisiTrace RS485 225 have, in reality, a shaft length of 215 mm. This ensures optimal rinsing in retractable armatures, such as Retractex.

ODO Cap L1: The L1 cap is designed for trace level measurements of dissolved oxygen in breweries, water de-aeration and power plants

Accessories



- **ODO Cap L1 Kit** Ref 10107102
- Calibration station Ref 243575

Cables see page ▶ 108 **Housings** see page **▶** 126

GIL

SENSORS



VisiTrace mA family



Specifications 0 to 2000 ppb (DO) Measuring range Measurement Principle Oxygen dependent luminescence quenching Response time t_{98%} < 20 s in gas; < 90 s in water Process temperature -20 to 140 °C, the sensor provides no DO reading above 85 °C Operating voltage 18 to 30 VDC Pressure range -1 to 12 bar (relative to ambient) Hygienic aspects Autoclavable, CIP, SIP Surface Quality $R_a < 0.4 \mu m (N5)$ Material Stainless steel 1.4435 O-ring **EPDM**

For more specifications see www.hamiltoncompany.com

The VisiTrace mA is designed to measure dissolved oxygen in the low ppb ranges in brewing applications, notably during filtration, and filling. In addition, the special designed ODO Cap L1 for breweries is stabilized against standard disinfectant solution with active chlorine and chlorine dioxide. This is powerful during measurements in breweries, which may not allow for calibration after every CIP.

With the transmitter integrated, the intelligent VisiTrace mA sensor provides more reliable measurements directly to your process control system via the 4-20 mA output. The also integrated Bluetooth 5 wireless interface may be used for monitoring, configuration and calibration, and saves time without compromising quality.



Did you know... that the VisiTrace mA is the only optical DO sensor that withstands chlorine and chlorine dioxide for a long time?

Benefits

- ► For measurements from 0 to 2000 ppb
- ► Stable against chlorine and chlorine dioxide
- ► Rapid start-up with no polarization
- ► Flow and CO₂ independent readings
- ► Robust design for high flow rates

Typical applications

- Brewerie
- Power Plante

Ordering Information

10068709								
•••••	Code	Interfac	се					
	1	mA/HAF	RT					
		Code	a-length	n (mm)				
		1	120					
		2	225*					
		3	325					
		4	425					
			Code	ODO Ca	ip .			
			1	L1				
				Code	Wetted Parts			
	+	+	+	1	EPDM			
10068709 –					← Order Code			

*The VisiTrace mA 225 have, in reality, a shaft length of 215 mm. This ensures optimal rinsing in retractable armatures, such as Retractex.

ODO Cap L1: The L1 cap is designed for trace level measurements of dissolved oxygen in breweries, water de-aeration and power plants.

Accessories



- ODO Cap L1 Kit Ref 10107102
- Calibration station Ref 243575

Cables see page ≥ 108

Housings see page ≥ 126

HART O







SENSORS

DO

VisiWater DO P



Specifications

Measuring range 0 to 40 ppm (DO)

Response time t_{98%} < 60 s at 25 °C, from air to nitrogen

Process temperature 0 to 60 °C

Pressure range -1 to 12 bar


Material Shaft: PVC-U
Cap: PPA

For more specifications see www.hamiltoncompany.com

The VisiWater DO P is an optical dissolved oxygen sensor designed for applications in water, wastewater, fish farming, lakes, and rivers. Its robust plastic shaft is ideal for these applications. The optical measurement technology ensures fast response time and minimum maintenance without polarization time. Like for all optical DO sensors the only spare part is the cap, which is easy and quickly replaceable.

The output signals 4-20 mA or Modbus can easily be integrated into process control systems (PCS). Calibration and configuration can be done via the PCS or ArcAir Desktop version with the help of the USB RS485 Modbus Converter.

Ordering Information



Benefits

- ➤ Simple and low maintenance
- ► Robust design
- ► Outdoor use incl. submersion

Typical applications

- Water and Wastewate
- Fish farming



Accessories



- **ODO Cap H20** Ref 243536
- Junction Box Ref 10076282

Cables see page ≥ 108

SENSORS **SENSORS**

OxyFerm FDA



The OxyFerm FDA is an electrochemical oxygen sensor suited for applications with high demands for hygiene, e.g. in pharmaceutical industry, in biotechnology and in food & beverage production. It is available with 12 mm or 25 mm (XL) shaft diameter.

The sensor is equipped with an FDA-approved membrane for use in hygienic processes. It withstands steam sterilization, autoclavation and CIP cleanings.

Benefits

- ► Sanitary Feature: The silicone membrane seals without a gap to steel membrane body (no additional o-ring)
- ► Little drift, fast response, short polarization time
- ► Replacing the cathode is possible and very simple to perform.

Typical applications





-670 mV

EPDM

For more specifications see www.hamiltoncompany.com

Ordering Information





Polarization voltage

O-ring





	a-length	T82	VP 6	Arc	MS
OxyFerm FDA	120	237450	237540	243100	237713
	160	237455	237541	243101	10069701
	225	237452	237542	243102	237715
	325	237453	237543	243103	10069700
	425	237454	237544	243104	_
OxyFerm XL	56	237175-OP	-	243140-OP	-
	125	237170	-	-	-
	262	237174	-	_	_
OxyFerm CIP	120	243289	-	-	-

With the XL option, the o-ring position can be optimally matched to the weld-in socket from 22 to 55mm. Please state the OP you need when ordering.

Accessories



- Membrane Kit FDA Ref 237140
- Membrane Kit CIP Ref 237126
- Membrane Kit Ref 237123
- Replacement Cathode OxyFerm Ref 237306
- Autoclavation Cap Oxyferm Ref 242000
- Polarization Module G Ref 237350
- Polarization Module T Ref 237370

Cables see page ≥ 108

Arc Accessories see page **▶** 117 **Housings** see page ▶ 126

93

92





• Oxylyte 30 mL Ref 237118

SENSORS



OxyGold B



Specifications Measuring range 8 ppb to 40 ppm (DO) Response time t98% < 60 s at 25 °C, from air to nitrogen 0 to 100 °C Process temperature Pressure range 0 to 12 bar (relative to ambient) CIP Hygienic aspects Oxylyte B Electrolyte Surface Quality $R_a < 0.4 \mu m (N5)$ Current in air at 25°C 180 to 500 nA Material Stainless steel 1.4435 Polarization voltage 0 mV **EPDM** O-ring

For more specifications see www.hamiltoncompany.com

The OxyGold B is an electrochemical oxygen sensor especially designed for applications which contain carbon dioxide like the production of beer, sparkling wine or soft drinks. The sensor is not affected by acidic gases.

Apart from the production of sparkling beverages, the OxyGold B can be used in all production processes where CO₂ might be an issue for electrochemical sensors.



Benefits

- ► No cross-sensitivity with CO₂
- Only very little flow required
- ► Pressure and CIP resistent
- ► Replacing the cathode is possible and very simple to perform.

Typical applications

- CO₂ recovery
- ▶ Water de-aeration



Ordering Information

	a-length	VP 6	Arc
OxyGold B	120	237180	not available
	225	237185	anymore*

*See VisiTrace sensor, page 92

Accessories



- OxyGold Membrane Kit Ref 237135
- Oxylyte B 30 mL Ref 237138
- Polarization Module B Ref 237360
- Replacement Cathode OxyGold B Ref 237437

Cables see page ≥ 108
Housings see page ≥ 126





SENSORS



OxyGold G



Specifications Measuring range 1 ppb to 40 ppm (DO) Response time t98% < 60 s at 25 °C, from air to nitrogen Process temperature 0 to 130 °C (Arc: analog 0 to 110 °C, digital 0 to 130 °C) Pressure range 0 to 12 bar (relative to ambient) Hygienic aspects Autoclavable, CIP, SIP Electrolyte Oxylyte G Surface Quality $R_a < 0.4 \mu m (N5)$ Current in air at 25°C 180 to 500 nA Stainless steel 1.4435 Material -670 mV Polarization voltage **EPDM** O-ring

For more specifications see www.hamiltoncompany.com

The OxyGold G is an electrochemical oxygen sensor designed for processes in which very small amounts of oxygen have to be traced, like in the pharmaceutical or microelectronics industry. It is also suitable for processes where high pressures are applied.

Benefits

- Trace level measuremen
- ➤ Suitable for use at high temperatures and high pressures during sterilization and CIP
- ► Little flow sensitivity
- ► Replacing the cathode is possible and very simple to perform.

Typical applications

- Boiler Feed Water
- Microelectronics



	a-length	VP 6	Arc
OxyGold G	120	237395	243110
	225	237396	243111

Accessories



- OxyGold Membrane Kit Ref 237135
- Oxylyte G 30 mL Ref 237139
- Polarization Module G Ref 237350
- Replacement Cathode OxyGold G Ref 237427

Cables see page ▶ 108

Arc Accessories see page **▶** 117

Housings see page ▶ 126







96 97

Ordering Information



Oxysens



Specifications	
Measuring range	40 ppb to 40 ppm (DO)
Response time t98%	< 60 s at 25 °C, from air to nitrogen
Process temperature	0 to 60 °C
Pressure range (relative to ambient)	0 to 4 bar
Electrolyte	Oxylyte
Surface Quality	R _a < 0.8 μm (N6)
Current in air at 25°C	40 to 80 nA
Material	Stainless steel 1.4435
Polarization voltage	-670 mV
O-ring	EPDM

For more specifications see www.hamiltoncompany.com

The Oxysens is an electrochemical oxygen sensor designed for applications in water, e.g. wastewater treatment, swimming pools or fish farms. It is easy to maintain, because the membrane and the electrolyte do not need to be replaced.

The response time of the Oxysens is fast, it is almost independent to flow and insensitive to soiling.

Benefits

- ➤ Maintenance-free DO sensor, no change of membrane or electrolyte
- ► Robust design
- ► Insensitive to soiling
- ► Short polarization and response times

Typical applications

- Water and Wastewater
- Fish farming



Ordering Information

	a-length	5 m fixed cable
Oxysens	120	237150
-		

Accessories



• Immersing Set Ref 237158

The Immersing Set sheaths and protects 120mm sensors such as Oxysens while immersed in streams or channels.

Housings see page ▶ 126





DURACAL pH BUFFERS

DURACAL pH BUFFERS

Buffer Solutions you can Trust

All calibration procedures assume that the labeled values of the calibration buffers are correct. But buffer values can change over time and so can your results. A complete range of patented buffer solutions provides pH stability up to 5 years, something never achieved before. The pH buffers 9.21 and 10.01 are even stable when exposed to air. High buffering capacity provides rapid, stable calibration. The growth of fungus and micro-organisms is prevented.



Traceability

An important issue for the production of Certified Reference Materials is to ensure traceability through an unbroken chain of comparisons to reference material of the highest metrological quality (Primary Reference Material) from NIST¹ and PTB². Unlike other manufacturers, where only top-down traceability is applied, Hamilton works with circular or closed-loop traceability, providing unique reliability of Hamilton DuraCal buffers.

Features

- ► Convenient 250 mL or 500 mL bottle with built-in calibration compartment
- ► Economical, only about 15 mL of buffer is used per calibration
- ➤ Certified pH value from a DAkkS laboratory accredited for pH measurement
- ► First class certificate with traceability to international standards
- ► Certificates available at www.hamiltoncompany.com
- Expiration date on the bottle
- ► Immune to microbial growth

Top-down traceability: At Hamilton, the pH value of DuraCal buffers is determined by comparison against two secondary reference buffer solutions from accredited suppliers of secondary reference materials. The solutions themselves are compared against primary reference solutions from PTB or NIST. The measurement uncertainties of every measurement comparison are known and documented.

Bottom-up traceability: To ensure the highest possible accuracy and full reliability of the pH value, a representative number of samples from every single production lot is verified by an external, independent and impartial DAkkS³ laboratory. The DuraCal samples are compared against secondary reference solutions from DAkkS and these are referenced themselves to primary reference solutions from PTB or NIST. At this stage, the traceability loop is closed. DAkkS provides Hamilton with a calibration certificate for every DuraCal production batch.

Certified reference material: Due to the complete traceability of the measurement procedure and the assignment of uncertainties to the particular testing steps, the buffers pH 4.01, 7.00, 9.21 and 10.01 are classified as "Certified Reference Material" (CRM).

pH Buffers

1.09 1.68 2.00 3.06 4.01 4.01 4.01 4.01 4.01 4.01 5.00 6.00 7.00 7.00 7.00 7.00 7.00 7.00 7	±0.02 ±0.02 ±0.02	60 60	Hamilton	500 mL	238271
2.00 3.06 4.01 4.01 4.01 4.01 4.01 4.01 5.00 6.00 7.00 7.00 7.00 7.00 7.00 7.00		60			
3.06 4.01 4.01 4.01 4.01 4.01 4.01 5.00 6.00 7.00 7.00 7.00 7.00 7.00 7.00	+0.02	00	Hamilton	500 mL	238272
4.01 4.01 4.01 4.01 4.01 4.01 5.00 6.00 7.00 7.00 7.00 7.00 7.00	±0.02	60	Hamilton	500 mL	238273
4.01 4.01 4.01 4.01 4.01 5.00 6.00 7.00 7.00 7.00 7.00 7.00	±0.02	60	Hamilton	500 mL	238274
4.01 4.01 4.01 4.01 5.00 6.00 7.00 7.00 7.00 7.00 7.00	±0.01/±0.02	24/60	DAkkS	250 mL	238317
4.01 4.01 4.01 5.00 6.00 7.00 7.00 7.00 7.00 7.00	±0.01/±0.02	24/60	DAkkS	500 mL	238217
4.01 4.01 5.00 6.00 7.00 7.00 7.00 7.00 7.00	±0.01/±0.02	24/60	DAkkS	3 x 500 mL	238917
4.01 5.00 6.00 7.00 7.00 7.00 7.00 7.00 7.00	±0.01/±0.02	24/60	DAkkS	5 L	238332
5.00 6.00 7.00 7.00 7.00 7.00 7.00	±0.01/±0.02	24/60	DAkkS	10 L	238194
6.00 7.00 7.00 7.00 7.00 7.00	±0.01/±0.02	24/60	DAkkS	1000 L	238895
7.00 7.00 7.00 7.00 7.00	±0.02	60	Hamilton	500 mL	238275
7.00 7.00 7.00 7.00	±0.02	60	Hamilton	500 mL	238276
7.00 7.00 7.00	±0.01/±0.02	24 / 60	DAkkS	250 mL	238318
7.00 7.00	±0.01/±0.02	24 / 60	DAkkS	500 mL	238218
7.00	±0.01/±0.02	24 / 60	DAkkS	3 x 500 mL	238918
•	±0.01/±0.02	24 / 60	DAkkS	5 L	238333
7.00	±0.01/±0.02	24 / 60	DAkkS	10 L	238188
	±0.01/±0.02	24 / 60	DAkkS	1000 L	238896
8.00	±0.02	60	Hamilton	500 mL	238277
9.21	±0.02	60	DAkkS	250 mL	238319
9.21	±0.02	60	DAkkS	500 mL	238219
9.21	±0.02	60	DAkkS	3 x 500 mL	238919
9.21	±0.02	60	DAkkS	10 L	238216
9.21	±0.02	60	DAkkS	1000 L	238897
10.01	±0.02	60	DAkkS	250 mL	238321
10.01	±0.02	60	DAkkS	500 mL	238223
10.01	±0.02	60	DAkkS	3 x 500 mL	238923
10.01	±0.02	60	DAkkS	10 L	238187
10.01	±0.02	60	DAkkS	1000 L	238898
11.00	±0.05	24	Hamilton	500 mL	238278
12.00	±0.05	24	Hamilton	500 mL	238279
4.01/7.00/9.21	±0.01/±0.02	24/60	DAkkS	500 mL, mixed	238922
4.01/7.00/10.01	±0.01/±0.02	24/60	DAkkS	500 mL, mixed	238924

Simple handling for professional results

Step 1 Open bottle



Step 2 Fill calibration compartment



Step 3 Calibrate electrode



Step 4 Empty calibration compartment

ORP Buffers

Value	Accuracy	Stability*	Certified By	Packaging Unit	Ref
271 mV	±5 mV	24	None	500 mL	238228
475 mV	±5 mV	24	None	250 mL	238322
475 mV	±5 mV	24	None	500 mL	238227

*In months after date of manufacturing

¹⁾ NIST: National Institute of Standards and Technology, Gaithersburg, MD, USA 2) PTB: Physikalisch Technische Bundesanstalt, Braunschweig, Germany

³⁾ DAkkS: Deutsche Akkreditierungsstelle GmbH (D-K-15186-01-00), Zentrum for Messen und Kalibrieren GmbH, Wolfen, Germany

CONDUCTIVITY STANDARDS

CONDUCTIVITY STANDARDS

Hamilton Conductivity Standards

Long-term stability and accuracy

For measurements in the low conductivity range stable and reliable calibration standards have been completely lacking up to now. Since a conductivity standard is not a buffer solution, the lower the value of the conductivity standard, the greater the effect of entry of CO2 or contamination. Hamilton is the first manufacturer to offer patented conductivity standards of 1.3 and 5 μ S/cm with a certified accuracy of $\pm 1\%$ and a lifetime of 1 and 3 years, respectively. The procedure for determining conductivity was developed in collaboration with DFM1. Many metrological institutes choose Hamilton standards because of their unprecedented stability and independent verification by PTB. During an interlaboratory test among prestigious European metrological institutes (PTB, DFM, DAkkS3) Hamilton standards were used as measurement solutions.



Hamilton is Different

Hamilton offers conductivity standards whose stability of $\pm 1\%$ is guaranteed over a lifetime of up to 3 years. They can be used repeatedly under the condition that the bottle is not left open for more than 1 hour in total.

A representative number of bottles from every batch are measured by DFM. Their value is recorded on the calibration certificate and on every bottle. DFM enjoys the highest prestige in Europe in the area of electrolytic conductivity and is equipped with an absolute measurement cell that was developed in collaboration with NIST, and is accredited by the Danish accreditation agency DANAK to a conductivity of 0.9 μS/cm. DFM and NIST⁴ have made comparisons of their measurement uncertainty and have confirmed in a series of scientific publications that the measurement accuracy is in each case the

same. Because no primary standards exist in the low conductivity range, measurements depend on absolute measurement cells which trace electrical conductivity back to the SI units: meter and volt. Testing of Hamilton standards is thus carried out on the most precise measurement apparatus in the world, and certified accordingly.




- 1) DFM: Danish Institute of Fundamental Metrology, Dänemark
- 2) PTB: Physikalisch-Technische Bundesanstalt, Braunschweig
- 3) DAkkS: Deutsche Akkreditierungsstelle4) NIST: National Institute of Standards and
- NIST: National Institute of Standards and Technology, Gaithersburg MD, USA

Unique advantages:

- Pemains stable for a minimum of 1 year for 1.3 μS/cm, and up to 3 years for all other values
- Certificate with calibration document from DFM (available at www.hamiltoncompany.com)
- Expiration date shown on every bottle
- ▶ Bottles are permitted to stay open for a total of 60 minutes

Stability of the Hamilton 5µS/cm Conductivity Standard over 36 months

Check measurement by PTB²







Value at 25°C	Accuracy	Stability*	Certificate From	Packaging Unit	Volume	Ref
1.3 μS/cm	±1%	12	DFM	Glass bottle	300 mL	238973
5 μS/cm	±1%	36	DFM	Glass bottle	300 mL	238926
15 μS/cm	±1%	36	DFM	Glass bottle	300 mL	238927
84 μS/cm	±1%	18	DFM	Calpack bottle	500 mL	238984
100 μS/cm	±1%	36	DFM	Glass bottle	300 mL	238934
147 μS/cm	±1%	18	DFM	Calpack bottle	500 mL	238985
706 μS/cm	±2%	36	Hamilton	Glass bottle	300 mL	238929
1413 µS/cm	±1%	36	DFM	Glass bottle	300 mL	238928
1413 μS/cm	±1%	18	DFM	Calpack bottle	500 mL	238986
12880 μS/cm	±1%	18	DFM	Calpack bottle	500 mL	238988
100 mS/cm	±1%	36	DFM	Glass bottle	300 mL	238935

*In months after date of manufacturing

ACCESSORIES

Oxygen Accessories



OxyFerm Membrane Kit

The OxyFerm Membrane Kit contains 3 membrane bodies, Oxylyte electrolyte, pipette, spare o-ring and a polishing strip.



Membrane Kit FDA

The Membrane Kit FDA is the kit for the OxyFerm FDA sensors and contains 3 FDA membrane bodies, Oxylyte electrolyte, pipette, spare o-ring and a polishing strip. The membrane body of the FDA membrane has a special rounded design to prevent accumulation of gas bubbles.



Membrane Kit CIP

The Membrane Kit CIP contains 3 membrane bodies that are especially designed to withstand CIP cleanings. Oxylyte electrolyte, pipette, spare o-ring and a polishing strip.

Ref	237126

OxyGold Membrane Kit

The OxyGold Membrane Kit contains 3 membrane bodies with the rounded design, pipette and a spare o-ring. Electrolyte must be ordered separately to match the sensor (see page 105).

Ref 237135

Polarization Module

The Polarization Module is to prepare replacement sensors so that they can be used immediately for measurements without connection to a transmitter. It polarizes the oxygen sensors and saves polarization time at the transmitter.

Polarization Module T OxyFerm / OxyFerm FDA / OxyFerm XL	Ref 237370
Polarization Module G OxyFerm VP / OxyGold G	Ref 237350
Polarization Module B OxyGold B	Ref 237360

Replacement Cathode OxyFerm	Ref 237306
Replacement Cathode OxyGold G	Ref 237427
Replacement Cathode OxyGold B	Ref 237437

Autoclavation Cap

The Autoclavation Cap is used to protect the OxyFerm T82 connector from moisture during autoclavation. It is important to keep connections dry and clean to ensure reliable measurements.

Autoclavation Cap OxyFerm Ref 242000

Electrolytes and Solutions



Electrolyte

3 M KCI	100 mL	238036
3 M KCI	500 mL	238936
Skylyte-CL	100 mL	242080
Protelyte	100 mL	238038
3 M KCI-LR	500 mL	238939
Skylyte	500 mL	238937
Electrolytes for Oxyger	Sensors	Ref
OxyGold Oxylyte G	30 mL	237139
OxyGold Oxylyte B	30 mL	237138
OxyFerm Oxylyte	30 mL	237118

Storage Solution

In order to to achieve long sensor life and faster electrode response times, it is recommended to store electrodes in our storage solution. It is an acid-buffered solution that ensures the regeneration of the electrode in addition to provide an optimized storage.

Storage Solution	500 mL	Ref 238931

ori Count Court ori Court ori Count Court ori Count Court ori Count Court ori C

Cleaning Solution Set

Depending on the type of application, the pH glass or diaphragm can get contaminated through various ingredients of the measuring solution. This is indicated by a slow response of the electrode, or even incorrect readings. To overcome these problems, Hamilton has developed a cleaning solution set. The intention is to have an overall cleaning of the pH glass as well as the diaphragm. The set is comprised of Cleaning Solution A, Cleaning solution B and a storage solution. To clean the electrode put it into each solution for 15 – 30 minutes, and your electrode will be ready for new measurements again.

Cleaning Solution Set Ref 238290


HAMILTON

Connectivity

An overview

Where and why, we need all these accessories

A quality measurement is nothing without a quality connection to your system. Whether a traditional analog connection or digitally via Modbus RS 485, we offer a broad range of connectivity options for you to choose from. The below diagram should help you navigate through the necessary requirements with ease.



Cables

A high quality measurement requires a high quality connection to the process control system. Hamilton cables ensure the best possible connection between your sensor and your process control system.

Sensor connector and relevant cables

SENSOR CONNECTION So what connector does my sensor have and what cable do I use? Below are a list of connectors available with Hamilton sensors.

CABLE CONNECTOR

The VP (VarioPin) is a common connector used throughout the Hamilton ensor product line. VP is abbreviation for "VarioPin". The VP designation often includes a number referring to the number of exposed.



K8 connectors are typically used on traditional pH / ORP sensors which lack temperature compensation. These connectors have a two pole design comprised of the center core and outer metallic threaded connection.



S7 and S8 connectors are typically found on traditional pH sensors with no temperature compensation. They are the same basic design however These connectors are recessed thus care must be taken to avoid moisture



Improved Electrical Properties

Indicator

Arrows

Robust **Design**

> **Hamilton** Logo






new

Introducing the **Hamilton made VP** connector

Now on all of our VP cables

Traditionally, VP connecters were every difficult to connect and disconnect. Our new connector was developed with special focus on the ease of connection.



Closing:

- Easy self alignment
- Closed position feedback

Opening:

- Tool less
- Low force



Cables for Traditional Sensors



or devices with open wiring terminals.


For sensors with standard (S7) connector. Controller side no connector (open end). Best suited for use with transmitters



For sensors with standard (S7) connector. Controller side

BNC connector. BNC connectors are commonly found on

Applikon biocontrollers and some older transmitters.





Length	Diameter	Ref
1 m	3 mm	355043
3 m	3 mm	355057
5 m	3 mm	355056





Diameter	Ref
3 mm	355045
3 mm	355059
	3 mm



For sensors with K8 connector. Controller side no connector (open end). Best suited for use with transmitters or devices with open wiring terminals.



Length	Diameter	Ref
1 m	5 mm	355153
3 m	5 mm	355154
5 m	5 mm	355155
10 m	5 mm	355156

For sensors with K8 connector. Controller side DIN connector. The DIN connector may be found on older Satorius biocontrollers and some laboratory pH meters.



Length	Diameter	Ref
1 m	5 mm	355157
2 m	5 mm	355158
3 m	5 mm	355159



For sensors with T82/D4 connector, e.g. OxyFerm. Controller side no connector (open end).



Length	Diameter	Ref
1 m	5 mm	355087
3 m	5 mm	355088
5 m	5 mm	355089
10 m	5 mm	355311

Cables for Intelligent Sensors

Connection for Industrial Processes e.g. Production (see page 13)

For sensors with T82/D4 connector, e.g. OxyFerm. Controller side Lemo connector.



Length	Diameter	Ref
1 m	5 mm	355160
2 m	5 mm	355161
3 m	5 mm	355162
5 m	5 mm	355163





Memosens

For sensors with Memosens connector. Controller side no connector (open end).





Length	Diameter	Ref
3 m	6.3 mm	355350
5 m	6.3 mm	355351
10 m	6.3 mm	355352



Compatible with:

• VisiFerm RS485-ECS family

VP8

- pH Arc family
- Conducell 4UxF family
- ORP Arc Sensors
- Conducell UPW Arc Sensors
- eDO Arc Sensor (e.g. OxyFerm FDA Arc)

Ref	Description	Interface
355263	1 m Data Cable VP8 / Open End	4-20 mA/Modbus
355264	3 m Data Cable VP8 / Open End	4-20 mA/Modbus
355265	5 m Data Cable VP8 / Open End	4-20 mA/Modbus
355266	10 m Data Cable VP8 / Open End	4-20 mA/Modbus
355267	15 m Data Cable VP8 / Open End	4-20 mA/Modbus
355268	20 m Data Cable VP8 / Open End	4-20 mA/Modbus
355217	1 m Cable VP8 / Open End	ECS mode*
355218	3 m Cable VP8 / Open End	ECS mode*
355219	5 m Cable VP8 / Open End	ECS mode*
355220	10 m Cable VP8 / Open End	ECS mode*
355221	15 m Cable VP8 / Open End	ECS mode*
355222	20 m Cable VP8 / Open End	ECS mode*
10109026	1m Data Cable (4 wire)	VP8/Open End new
10109251	2m Data Cable (4 wire)	VP8/Open End new
10109250	3m Data Cable (4 wire)	VP8/Open End new



VP6

For sensors with VP6 connector. VP6 single coaxial cable. Controller side no connector (open end).



Length	Diameter	Ref
1 m	7,5 mm	355108
2 m	7,5 mm	355187
3 m	7,5 mm	355109
5 m	7,5 mm	355110
10 m	7,5 mm	355111
20 m	7,5 mm	355112



Compatible with:

• All Arc Sensors



Compatible with:

All Arc Sensors

Ref	Description
10070910	1 m Data Cable VP8 / M12-8 Pole (male)
10071905	3 m Data Cable VP8 / M12-8 Pole (male)
10067844	5 m Data Cable VP8 / M12-8 Pole (male)
10067846	10 m Data Cable VP8 / M12-8 Pole (male)

Ref	Description
10108609	1m Cable VP8 (F) / VP8 (F)
10108610	2m Cable VP8 (F) / VP8 (F)
10108611	3m Cable VP8 (F) / VP8 (F)

^{*} VisiFerm DO family only

CABLES

Power Cables for Bio Controllers


Connection for old Bio Controllers or Transmitters in R&D (see page 15)

If you want to gain the benefits our Arc Intelligent sensors can give you but need to stick with an analog sensor connection with your transmitter or PCS, the following cables can assit in giving you this backwards capability.



M12 8-Pole

Description
3 m Cable M12-8 Pole / Open End
5 m Cable M12-8 Pole / Open End
10 m Cable M12-8 Pole / Open End



Ref	Description
	· · · · · · · · · · · · · · · · · · ·
355298	1 m Power Cable VP8 / AMP
355258	4 m Power Cable VP8 / Binder
355297	1 m Power Cable VP8 / BNC
355296	3 m Power Cable VP8 / BNC
355245	2.5 m Power Cable VP8 / Lemo
355194	1 m Power Cable VP8 / Open End



M12 4-Pole



Compatible with:

- VisiFerm mA family
- VisiTrace mA familiy

M12-4	Pole	Ex Protection	230 V / 110 V AC
		HAMILTON	

Ref	Description
355283	3 m Cable M12-4 Pole / Open End
355284	5 m Cable M12-4 Pole / Open End
355285	10 m Cable M12-4 Pole / Open End

Ref	Description
355288	3 m Power Cable M12-4 Pole

VP8 Arc ECS Adapter pH/ORP BNC / Open End Compatible with: 230 V / 110 V AC **Power Cable** pH Arc family • ORP Arc family

For retrofit of existing polarographic DO sensor installations with VisiFerm RS485-ECS sensors.



Compatible with:

• VisiFerm RS485-ECS family

Ref	Description
243168-XX	Arc ECS Adapter pH/ORP BNC
243169-XX	Arc ECS Adapter pH/ORP Open End

The code XX in the product number defines the type of electrical

01 – Power cord EU / 02 – Power cord CH / 03 – Power cord US 04 - Power cord UK / 05 - Power cord AU/NZ

Ref	Description
242413-XX	VisiFerm T82/D4-Power Adapter

The code XX in the product number defines the type of electrical


01 – Power cord EU / 02 – Power cord CH / 03 – Power cord US 04 – Power cord UK / 05 – Power cord AU/NZ

Compatible with:

- VisiFerm mA family
- VisiTrace mA familiy

CABLES ARC ACCESSORIES


The new Power Cable M12-8 Pole / open end is designed for use with the Arc Wi 2G Adapter BT (Ref 243470) to facilitate an "active" 4-20 mA signal.



Ref	Description
10143091	1m Power Cable M12-8 Pole / open end / power plug new
10143092	3m Power Cable M12-8 Pole / open end / power plug new

Cables for connection to Arc Sensors

For connecting Arc sensors to ArcAir software



Ref	Description
243490-01	2 m Arc USB Power Cable VP8
243490-02	2 m Arc USB Power Cable M12-8 Pole
242176	2 m Arc Cable VP8 / M8
355339	2 m Service Cable M12-8 Pole / M8
355289	2 m Service Cable M12-4 Pole / M8

*For connection with the Arc USB power cable or Arc Modbus OPC Converter

Arc Accessories

Arc Wireless Converter BT





Designed for wireless communication between ArcAir and Arc sensors via Computer. ArcAir Advanced license included.

Ref 242333

USB RS485 Modbus Converter



Designed for wired communication between ArcAir and Visiwater DO fix cable sensor.

Ref 242411

Arc Wi 1G Adapter BT Arc Wi 2G Adapter BT



These Adapters are expanding the functionality of Arc sensors by providing wireless communication for local monitoring all analog and digital signals, in parallel to robust 4-20 mA signal, and simple sensor connection to the PCS.

Description		
Arc Wi 1G Adapter BT		
Arc Wi 2G Adapter BT		

Arc View Mobile



This mobile device empowers the operator to monitor measurement values, calibrate Arc sensors and configure various parameters with a unified user interface for pH, DO, Conductivity and ORP. The Arc View Mobile device is based on the Samsung Galaxy Tab Active tablet and comes pre-configured with the ArcAir application, app blocker application and power supply.

Description
Arc View Mobile Basic
Arc View Mobile Advanced

Digital Converters

Arc Modbus OPC Converter



The Arc Modbus OPC converter converts the Modbus protocol of Arc sensors into OPC UA protocol. All Arc parameters are supported.

10089359

Ref

Modbus Profinet Converter



The Arc Modbus Profinet Converter ist a 24VDC DIN rail mounted device that can power up to four Arc sensors and converts their Modbus output to Profinet.

Ref

10116586

Modbus Profibus Converter



The Arc Modbus Profibus Converter is a 24VDC DIN rail mounted device that powers up to four Arc sensors and converts their Modbus output to Profibus DP.

Ref

243555

Hamilton Customized Products

Customized products for our customers' special needs

The adaptation of standard products to customer's special needs is the main focus of our application engineering team. Customizing can include modifications to length, insertion depth, process adaptation of the sensor or changing the housing to a different material. Many more adaptions are possible.



HAMILTON' CUSTOMIZED PRODUCT

Need a custom housing or sensor? The Hamilton Customized Product Team is happy to help design products for your specific application. Give us a call to learn more.



Transmitter H100

рН

The H100 is a transmitter for universal use in the chemical industry, power stations, biotechnology, food processing and pharmaceutical industries as well as in water/wastewater treatment. Icons guide the operator and show the sensor status.

Sensor failures are detected, shown on the display and an alarm is set. Calibration can be done manually or by selecting standard calibration media. After each calibration the sensor data will be shown and evaluated. The H100 is easy to handle and can be mounted on the wall as well as on a panel.



User friendly, robust and reliable



Easy to install, operate and calibrate

- Large terminal compartment and pre-assembled rear unit for easy installation.
- ▶ The large display and intuitive menu structure ensure straightforward navigation.
- lcons supply operating messages and signal unusual states.
- Simple calibration with automatic buffer recognition.



Robust design

- Optional protective hood for additional protection against weather exposure and mechanical damage.
- Wall, post/pipe, or panel mounting possible with optional panel- or pipe-mount kit.



Reliable instrument for process applications

- The sensor status and potential defects are continuously monitored for real time display of error or alarm.
- Asymmetry potential, slope and response time are evaluated during calibration through the sensor lifetime for preventive maintenance indication.
- The integrated calibration timer automatically indicates when calibration is required.

Transmitter H100 pH

Measured variable	pH, mV and temperature
Measuring range (pH / OPR)	-1500 to +1500 mV
Display range pH	-2.00 to 16.00
Measuring error	<0.02 pH, <1 mV
Temperature input	Pt 100, Pt 1000, NTC 30 kOhm
Temperature measuring range	-20.0 to +150 °C
Temperature resolution	0.1 °C
Calibration	1 point, 2 point and product calibration
Power supply	24 to 230 V AC/DC
Display	LC display, 7-segment with icons
Ambient temperature	-20 to 55 °C
Relative humidity	80 % at temperatures up to 55 °C
Ingress protection	IP 65, NEMA 4X
Alarm contact	Yes
Hold mode	Yes
FM	NI Class I, Div 2 Group A, B, C, D, T4 Ta = 55 °C; Type 2 NI Class I, Zone 2 Group IIC, T4 Ta = 55 °C; Type 2

Ordering Information

Туре	Ref
H100 pH	243080-01




Accessories

- Pipe-mount kit Ref 243082
- Panel-mount kit Ref 243083
- Protective hood Ref 243084

Mounting plan

all dimensions in mm



- 1 Cable gland (3x)
- 2 Knockouts for cable glands or ½" conduit (conduits not incl.)
- 3 Knockout for pipe mounting (4x)
- 4 Knockout for wall mounting (2x)





Transmitter H100 Cond

Measured variable	Conductivity, resistivity, concentration, salinity, temperature
Measuring range conductivity	0 to 999.9 mS/cm
Effective range conductivity	0.2 μS x c to 1000 mS x c
Measuring range resistivty	0.00 to 99.99 MΩ x cm
Measuring range concentration	0.00 to 9.99 % by wt
Measuring range salinity	0.0 to 45 ‰ (0 to 35 °C)
Measuring error	< 1 % meas. val. + 0.4 µS x c
Temperature input	Pt 100, Pt 1000, NTC 30 kOhm
Temperature measuring range	Pt 100/Pt 1000: -20.0 to +200 °C NTC 30 kOhm: -20.0 to +150 °C
Temperature resolution	0.1 °C
Power supply	24 to 230 V AC/DC
Display	LC display, 7-segment with icons
Ambient temperature	-20 to 55 °C
Relative humidity	80 % at temperatures up to 55 °C
Ingress protection	IP 65, NEMA 4X
Alarm contact	Yes
Hold mode	Yes
FM	NI Class I, Div 2 Group A, B, C, D, T4 Ta = 55 °C; Type 2 NI Class I, Zone 2 Group IIC, T4 Ta = 55 °C; Type 2

Ordering Information

Туре	Ref
H100 Cond	243080-02



Accessories

- Pipe-mount kit Ref 243082
- Panel-mount kit Ref 243083
- Protective hood Ref 243084

Transmitter H100 DO

Measured variable	DO saturation, DO concentration
Measuring current	-2 to 1800 nA
O ₂ resolution	0.05 nA
O ₂ saturation	0 to 200 %
O ₂ concentration	0.00 to 20.00 mg/l / 0.00 to 20.00 ppm
Polarization voltage	0 to 1000 mV (User-defined)
Salinity correction	00.00 to 45.00 g/kg (User-defined)
Measuring error	< 0.5 % meas. val. + 0.5 %
Temperature input	NTC 22 kOhm, NTC 30 kOhm
Temperature measuring range	-20.0 to +150 °C
Temperature resolution	0.1 °C
Power supply	24 to 230 V AC/DC
Display	LC display, 7-segment with icons
Ambient temperature	-20 to 55 °C
Relative humidity	80 % at temperatures up to 55 °C
Ingress protection	IP 65, NEMA 4X
Alarm contact	Yes
Hold mode	Yes
FM	NI Class I, Div 2 Group A, B, C, D, T4 Ta = 55 °C; Type NI Class I, Zone 2 Group IIC, T4 Ta = 55 °C; Type 2

Ordering Information

Туре	Ref	
H100 DO	243080-03	




Accessories

- Pipe-mount kit Ref 243082
- Panel-mount kit Ref 243083
- Protective hood Ref 243084


Mounting plan


all dimensions in mm



- 1 Cable gland (3x)
- ½" conduit (conduits not incl.)

Mounting plan





- 1 Cable gland (3x)
- 2 Knockouts for cable glands or ½" conduit (conduits not incl.)
- 3 Knockout for pipe mounting (4x)
- 4 Knockout for wall mounting (2x)

2 Knockouts for cable glands or

- 3 Knockout for pipe mounting (4x)
- 4 Knockout for wall mounting (2x)

Transmitter H220X









Hamilton H220X Transmitters combine ease of use and reliability. They are available in different configurations: Analog pH / ORP, Conductivity and inductive Conductivity as well as Memosens® pH and Oxygen.

It has been designed for universal process application including use in pharmaceutical, chemical, food & beverage industries as well as water / waste water treatment. The self-explaining user interface ensures comfortable and intuitive handling. Hamilton H220X transmitters provide continuous sensor monitoring and preventive maintenance indication for maximal reliability. The Memosens® Technology allows plug & play with pre-calibrated Memosens® sensors. Predictive maintenance system detects when a sensor has to be cleaned, calibrated or replaced.



Perfectly designed for hazardous areas and the Memosens® technology



Easy to install, operate and calibrate

- The large display and intuitive menu structure ensure straightforward navigation
- ► Simple calibration with automatic buffer recognition
- Memosens® sensors can be connected for even more simple handling

Robust design

- Suitable for Explosions protected areas (Ex II (1) 2G Ex ib [ia Ga] IIX T6/T4 Gb)
- ► Wall, post/pipe, or panel mounting possible
- Transmitter suitable for pollution degree 3




Reliable instrument for process applications

- Sensor status and potential defects are continuously monitored; errors and alarms are displayed in real time
- Asymmetry potential, slope and response time are evaluated during calibration through the sensor lifetime for preventive maintenance indication
- User-guided commissioning, graphic display and plain text guidance for maximum operating safety


 pH / ORP analog 	
• pH / ORP Memosens	
Conductive Conductivity	analog
 Inductive Conductivity a 	nalog
• eDO Memosens	
More info about measuring on the Hamilton website.	ranges, temperature ranges, input and output signals can be found

243081				
	Code	Sensor	Module	
	1	Conductivity, Conductive Sensor		
	2	Conductivity, Inductive Sensor		
	3	Digital, Memosens pH, ORP		
	4	Digital, N	Memosens eDO	
	5	pH or Of	RP (analog)	
		Code	Software	
		1	Standard Version	
	+	2	Advanced Version	
243081 -			← Order Code	

Mounting plan



220



all dimensions in mm

Front —

— Bottom ——



Different processes have different requirements for sensors to provide an accurate and reliable measurement. Being in contact with the media is the most important one. In order to meet the different requirements, Hamilton has developed various kinds of housings and armatures: static, retractable, pressurizable, pneumatic, manual, weld-in and hygienic sockets.

No matter what type of housing is needed for a pipe or a vessel, on the following pages the right one for each application can be found.

FlexiFit



The FlexiFit housings are designed for 120 mm sensors in different kinds of industries. A variety of process connections ensure the usability in the chemical industry as well as in hygienic processes. All FlexiFit have EPDM o-rings and the electropolished surface quality (Ra < 0.4 μ m) quality is shown on a certificate. They are suitable for autoclavation, CIP and SIP procedures.

Benefits

- ► Easy installation and handling
- ➤ Various o-ring positions available
- ► Hygienic design

Ordering Information

Process Connection	Ref
G 11/4	237331-OP
G 1¼	237380-OP
TC 1.5"	237341
Varivent®	237344
Varivent®	237345
	G 11/4 G 11/4 TC 1.5" Varivent®





U = Unprotected TC = Triclamp


Specifications Wetted parts Stainless Steel 1.4435 O-ring material **EPDM** O-ring position 22 to 55 mm (G 11/4) Pressure range (relative to ambient) 0 to 6 bar -10 to 140 °C Temperature range Sensor thread PG 13.5 Sensor a-length 120 mm Surface finish $R_a < 0.4 \mu m$ (N5 electropolished)

For more specifications see www.hamiltoncompany.com


Other designs and materials available on request

Dimensional drawing / FlexiFit

all dimensions in mm



FlexiFit VV-0



FlexiFit Bio

Accessories



- Service Kit FlexiFit Bio Ref 237366
- Service Kit FFKM Ref 237319
- Service Kit FKM Ref 237219
- Service Kit RetractoFit PEEK Ref 237388

Safety Socket see page ▶ 152



The RetractoFit is a retractable armature designed for 225 mm sensors in industrial applications. It allows the operator to mount and dismount sensors while the process is running. Safe sensor handling during process is guaranteed because insertion into the vessel without a sensor is impossible so is removal while in the measuring position. It is easy to use and maintain: only one press on the red button is needed to move the sensor into or out of the process. All o-rings can easily be replaced by the operator without special tools. The RetractoFit is available in different versions.

When the housing with an Arc sensor, VisiFerm mA, VisiTrace mA and protective sleeve the aperture (hole) in the protective sleeve must be enlarged or the housing has to be used without the protective sleeve. Wireless adapters on top of Arc sensors can only be used without the protective sleeve.





Benefits

- ► Integral safety mechanism
- ➤ Sensor can be withdrawn from the process for cleaning, calibration or replacement
- ► Easy maintenance
- ➤ 3.1 material certificate included

Ordering Information


Туре	Process Connection	Ref
RetractoFit	G 11/4	237240
RetractoFit PEEK 25	G 1¼	237490

Specifications Wetted parts RetractoFit: Stainless Steel 1.4571 RetractoFit PEEK: PEEK (FDA approved) O-ring material O-ring position RetractoFit: 22.5 mm RetractoFit PEEK: 25 mm Pressure range (relative to ambient) 0 to 6 bar -10 to 130 °C Temperature range PG 13.5 Sensor thread Sensor a-length 225 mm Surface finish RetractoFit: Ra < 0.4 µm (N5 electropolished)


For more specifications see www.hamiltoncompany.com

Dimensional drawings / RetractoFit

all dimensions in mm



Maintenance position



Measuring position

Accessories



- Service Kit RetractoFit Ref 237239
- FFKM Kit RetractoFit Ref 237339

Safety Socket see page ▶ 152

RetractoFit Bio



The RetractoFit Bio is a retractable housing designed for 225 mm sensors in hygienic applications in the biotechnology, food & beverage and the pharmaceutical industry. It allows the operator to mount and dismount sensors while the process is running. Safe sensor handling during the process is guaranteed because insertion into a vessel without sensor is impossible so is removal while in the measuring position. It is easy to use and maintain: only one press on the red button is needed to move the sensor into or out of the process. All o-rings can be easily be replaced by the operator without special tools.




Did you know... that the RetractoFit Bio has a special rinsing chamber with angled connections for cleaning solutions and special inlet construction guarantees an entire cleaning of the chamber through a swirl effect



Benefits

- ► Integral safety mechanism
- ➤ Sensor can be withdrawn from the process for cleaning, calibration or replacement
- Special hygienic design of cleaning chamber
- ► Easy maintenance




Specifications Wetted parts Stainless Steel 1.4435 O-ring material **EPDM** 22 mm and 55 mm O-ring position Pressure range (relative to ambient) 0 to 6 bar -10 to 140 °C Temperature range Sensor thread PG 13.5 Sensor a-length 225 mm R_a < 0.4 µm (N5 electropolished) Surface finish

For more specifications see www.hamiltoncompany.com

Dimensional drawings / RetractoFit Bio 25

all dimensions in mm

Maintenance position



Measuring position

Ordering Information

Туре	Process Connection	Ref
RetractoFit Bio 25	G 11/4	237480
RetractoFit Bio 55	G 1¼	237440

Accessories



• FDA Service Kit Ref 237338

Safety Socket see page ▶ 152








The retractable pneumatic or manual housing Retractex B was designed for sanitary applications in biotechnology, food & beverage and pharmaceutical industry. The compact design with a stroke of only 36 mm keeps wear on seals to a minimum and creates excellent reliability – day and night, all year long. It can be cleaned easily and thoroughly in place, including the space between the socket and rinsing chamber. The Retractex B with its patented HyCIP cleaning principle offers the best available cleaning efficiency for Ingold sockets (G 11/4").

It is designed for 12 mm sensors and is equipped with several safety features (e.g. no sensor – no insertion, window to check seals for leakage etc.) to provide operater safety. It is available with various process connections that can be used with all vessels used in these branches.



How does the HyCIP process connection work?

In cleaning position, the sensor can be cleaned and sterilzed together with all wetted seals. In the HyCIP connection the cleaning solution is directed between armature and socket up to the process seal so the most remote parts of the chamber are rinsed. Thus HyCIP housings are unmatched for their cleaning performance of the sensor and of all relevant seals.

Benefits

- ► Extremely compact design
- ► Integrated safety concept- no sensor no insertion
- ➤ Very low maintenance
- ► Sterile safety and unique cleaning efficiency with HyCIP



For more specifications see www.hamiltoncompany.com

Ordering Information

243240	Retractex B (pneumatic)										
243275	Retractex B M (manual)										
	Code	Materia	al (wetted	parts)							
	1	Stainles	s Steel 1.4	1404							
	0	special									
		Code	Sealing	g Material	l (wetted s	ealings)					
		1	EPDM/FDA USP class VI								
		2	FKM								
		0	special								
			Code	Sensor	r						
			1	225 mn	n PG13,5						
			0	special							
				Code	Process	s Connect	tion				
				1	Ingold (0	3 1¼") o-Ri	ng Position 28 mm				
				2	Varivent N DN 40-125						
				3		Ø 50,5 mm)					
				4	TriClamp 2" (OD Ø 64 mm)						
		5	6 DIN 11851 DN50								
	7				, , , , , , , , , , , , , , , , , , ,						
				8		HyCIP for Ingold (G 11/4") o-Ring Position 50 mm					
				9	HyCIP for Ingold (G 11/4") o-Ring Position 55 mm						
				0	special						
					Code		g Connection				
					1		read female				
					2		ead female				
					3	1/4" NPT					
					9		o ¾" Ø 4 mm o ¾" Ø 10,3 mm (Sartorius)				
							7 74 W 10,3 IIIII (Sartorius)				
					0	special	Position switch				
						Code 1	pneumatic / without for manual				
						0	special				
					<u> </u>		← Order Code				

Retractex BC Steel




The retractable pneumatic or manual housing Retractex BC is designed for applications in the chemical industry. The compact design with a stroke of only 36 mm keeps wear on seals to a minimum and creates to excellent reliability – day and night, all year long. It can be cleaned easily and thoroughly in place. It is designed for 12 mm sensors and is equipped with several safety features (e.g. no sensor – no insertion, window to check seals for leakage etc.) to provide operator safety. The Retractex BC comes with a G 11/4" process connection and is available with two different o-ring positions.

Cleaning of the Retractex BC?

In cleaning position, the sensor can be cleaned while the process is running. The advantage of the insertion tube is the short way for insertion. A PTFE scraper with o-ring guarantees that dirt stays outside of the armature and does not harm the o-ring.

Benefits

- Extremely compact design
- ▶ Integrated safety concept- no sensor no insertion
- Very low maintenance




For more specifications see www.hamiltoncompany.com

Ordering Information

237730	Retrac	Retractex BC Steel (pneumatic)											
237735		tex BC Ste											
	Code	Materia	ıl (wetted	parts)									
	1	Stainles	s Steel 1.4	404									
	2	Stainles	s Steel 2.4	602									
	0	special											
		Code Sealing Material (wetted sealings)											
		1	1 EPDM/FDA USP VI										
		2	FKM										
		3	FFKM										
		0 special											
			Code	Sensor									
			1	1 225 mm PG13,5									
			0	special									
				Code	Process Connection								
				1	Ingold (G 11/4") o-Ring Position 28 mm								
				2	Ingold (G 11/4") o-Ring Position 50 mm								
				0	special								
					Code		g Connection						
					1		ead female						
					2		ead female						
					3	1/4" NPT female							
					0	special							
						Code	Position switch						
						1	pneumatic						
		_			_	0	special						
Ref							← Order Code						

Retractex BC Plastic




The retractable pneumatic or manual housing Retractex BC was designed for applications in the chemical industry. The compact design with a stroke of only 36 mm keeps wear on seals to a minimum and creates excellent reliability – day and night, all year long. It can be cleaned easily and thoroughly in place. It is designed for 12 mm sensors and is equipped with several safety features (e.g. no sensor – no insertion, window to check seals for leakage etc.) to provide operator safety. The Retractex BC comes with a G 1½" process connection and is available with two different o-ring positions.

Cleaning of the Retractex BC?

In cleaning position, the sensor can be cleaned while the process is running. The advantage of the insertion tube is the short way for insertion. A PTFE scraper with o-ring guarantees that dirt stays outside of the armature and does not harm the o-ring.

Benefits

- ► Extremely compact design
- ► Integrated safety concept- no sensor no insertion
- ► Very low maintenance
- ► Easy installation of the pneumatic armature with color coded connectors
- ► Choice of 3 different plastics



For more specifications see www.hamiltoncompany.com

Ordering Information

237740	Retract	tex BC Plastic (pneumatic)									
237745		ex BC Pla									
	Code	Materia	I (wetted	l parts)							
	1	PP									
	2	PVDF / S	Stainless	Steel 2.460	2						
	3	PEEK									
	0	special									
		Code Sealing Material (wetted sealings)									
		1	EPDM/	FDA USP V	l						
		2	FKM								
		3	FFKM								
		0	special	special							
			Code	Code Sensor							
			1 225 mm PG13,5								
			0	special							
				Code							
				1	Ingold (G 11/4") o-Ring Position 25 mm						
				0	special						
					Code		ng Connection				
					1		read female				
					2		read female				
					3	1/4" NPT	female				
					0	special					
						Code	Position switch				
					1	pneumatic / without for manual					
	_	_	_			0	special				
Ref							← Order Code				


Retractex C Steel



The retractable pneumatic or manual housing Retractex C was designed for applications in the chemical industry. The compact design with a stroke of only 36 mm keeps wear on seals to a minimum and creates to excellent reliability – day and night, all year long. It can be cleaned easily and thoroughly in place. It is designed for 12 mm sensors and is equipped with several safety features (e.g. no sensor – no insertion, window to check seals for leakage etc.) to provide operator safety. It is available with various process connections that can be used with all vessels used in this branch.


Cleaning of the Retractex C?

In cleaning position, the sensor can be cleaned while the process is running. The advantage of the insertion tube is the short way for insertion. A PTFE scraper with o-ring guarantees that dirt stays outside of the armature and does not harm the o-ring.



Benefits

- ► Extremely compact design
- ▶ Integrated safety concept- no sensor no insertion
- ➤ Very low maintenance
- ► Easy installation of the pneumatic armature with color coded connectors



Stainless Steel 1.4404 or 2.4602

Ordering Information

Specifications
Wetted parts

243200 243255	Retractex C Steel (pneumatic) Retractex C Steel M (manual)											
	Code		l (wettec									
	1		s Steel 1.4									
	2		s Steel 2.4									
	0	special										
		Code	Sealin	g Material	(wetted s	ealings)						
		1	EPDM /									
		2	FKM									
		3	FFKM									
		0	special									
			Code	Sensor								
			1	225 mm	PG13,5	PG13,5						
			0									
				Code	Process Connection							
				1	Flange DN32 PN16							
				2	Flange DN40 PN16							
				3	Flange DN50 PN16							
				4	Flange A	ANSI 1¼" 1	50lbs					
				5	Flange ANSI 11/2" 150lbs							
				6		Flange ANSI 2" 150lbs						
				7 NPT M 1¼"								
				8	TriClam	TriClamp 2"						
				0 special								
					Code		g Connection					
					1		read female					
					2		ead female					
					3	1/4" NPT	female					
					0	special						
						Code	Position switch					
						1	pneumatic / without for manual					
	+	+	+	+	+	0	special					
Ref							← Order Code					



HOUSINGS HOUSINGS


Retractex C Plastic



The retractable pneumatic or manual housing Retractex C was designed for applications in the chemical industry. The compact design with a stroke of only 36 mm keeps wear on seals to a minimum and creates excellent reliability – day and night, all year long. It can be cleaned easily and thoroughly in place. It is designed for 12 mm sensors and is equipped with several safety features (e.g. no sensor – no insertion, window to check seals for leakage etc.) to provide operator safety. It is available with various process connections that can be used with all vessels used in this branch.


Cleaning of the Retractex C?

In cleaning position, the sensor can be cleaned while the process is running. The advantage of the insertion tube is the short way for insertion. A PTFE scraper with o-ring guarantees that dirt stays outside of the armature and does not harm the o-ring.



Benefits

- ► Extremely compact design
- ▶ Integrated safety concept- no sensor no insertion
- ➤ Very low maintenance
- ► Easy installation of the pneumatic armature with color coded connectors
- ► Choice of 3 different plastics



Ordering Information

Specifications Wetted parts

O-ring material

Temperature range

Sensor a-length

Surface finish

Pressure range (relative to ambient)

For more specifications see www.hamiltoncompany.com

243220	Retract	ex C Plastic (pneumatic)								
243265		ex C Plas								
	Code	Materia	I (wetted	parts)						
	1	PP								
	2	PVDF / Stainless Steel 2.4602 PEEK special								
	3									
	0									
		Code	Code Sealing Material (wetted sealings)							
		1	EPDM /	FDA USP	class VI					
		2	FKM							
		3	FFKM							
		0	special							
			Code	Sensor						
			1	225 mm						
			0	special						
				Code		s Connect	tion			
				1	Flange D					
				2	Flange ANSI 2" NPT M 1¼"					
				3						
				0	special	01 :	0 "			
					Code		ng Connection read female			
					2		read female			
					3	1/4" NPT				
					0	special				
						Code	Position switch			
						1	pneumatic / without for manual			
	1	1	 	1	<u> </u>	0	special			
Ref							← Order Code			

PVDF or PEEK or PP

-10 to 140 °C

 $R_a < 0.8 \ \mu m \ (N6)$

225 mm

EPDM or FKM or FFKM

0 to 16 bar (120 °C), 10 bar (140 °C)



Retractex C Steel LT




The retractable pneumatic or manual housing Retractex C was designed for applications in the chemical industry. The compact design with a stroke of only 36 mm with an insertion depth up to 207 mm keeps wear on seals to a minimum and creates excellent reliability – day and night, all year long. It can be cleaned easily and thoroughly in place. It is designed for 12 mm sensors and is equipped with several safety features (e.g. no sensor – no insertion, window to check seals for leakage etc.) to provide operator safety. It is available with various process connections that can be used with all vessels used in this branch.

Cleaning of the Retractex C?

In cleaning position, the sensor can be cleaned while the process is running. The advantage of the insertion tube is the short way for insertion. A PTFE scraper with o-ring guarantees that dirt stays outside of the armature and does not harm the o-ring.





Benefits

- ► Extremely compact design (only 36 mm trave of inertion tube with an insertion depth of 207 mm)
- ▶ Integrated safety concept- no sensor no insertion
- ➤ Very low maintenance
- ► Easy installation of the pneumatic armature with color coded connectors

SpecificationsWetted partsStainless steel 1.4404 or 2.4602O-ring materialEPDM or FKM or FFKMPressure range (relative to ambient)0 to 16 bar (120 °C), 10 bar (140 °C)Temperature range-10 to 140 °CSensor a-length325 mmSurface finishRa < 0.8 μm (N6)</td>

 $For \ more \ specifications \ see \ www.hamiltoncompany.com$

Ordering Information


243210	Retract	tex C Stee	LT (pne	umatic)					
243260		tex C Stee							
	Code	Materia	l (wetted	parts)					
	1	Stainless	Steel 1.4	404					
	2	Stainless	Steel 2.4	602					
	0	special							
		Code	Sealing	Material	(wetted s	ealings)			
		1	EPDM /	FDA USP	class VI				
		2	FKM						
		3	FFKM						
		0	special	•					
			Code	325mm PG13,5					
			1						
			0	special					
				Code		Connect	tion		
				1	Flange D				
				2	Flange D				
				3	Flange A				
				4	Flange A	ANSI 2"			
				0	special				
					Code		g Connection		
					1		read female		
					3	14" NPT 1	ead female		
					0		IEIIIdIE		
						special	Position switch		
						1	pneumatic / without for manual		
						0	special		
						J 0	opolici.		
Ref							← Order Code		

Retractex C
Plastic LT

The retractable pneumatic or manual housing Retractex C was designed for applications in the chemical industry. The compact design with a stroke of only 36 mm with an insertion depth up to 207 mm keeps wear on seals to a minimum and creates excellent reliability – day and night, all year long. It can be cleaned easily and thoroughly in place. It is designed for 12 mm sensors and is equipped with several safety features (e.g. no sensor – no insertion, window to check seals for leakage etc.) to provide operator safety. It is available with various process connections that can be used with all vessels used in this branch.

Cleaning of the Retractex C?

In cleaning position, the sensor can be cleaned while the process is running. The advantage of the insertion tube is the short way for insertion. A PTFE scraper with o-ring guarantees that dirt stays outside of the armature and does not harm the o-ring.



Benefits

- Extremely compact design (only 36 mm travel of insertion tube with an insertion depth of 207 mm)
- ► Integrated safety concept- no sensor no insertion
- ► Very low maintenance
- ► Easy installation of the pneumatic armature with color coded connectors

Wetted parts	PVDF or PEEK
O-ring material	EPDM or FKM or FFKM
Pressure range (relative to ambient)	0 to 16 bar (120 °C), 10 bar (140 °C)
Temperature range	-10 to 140 °C
Sensor a-length	325 mm
Surface finish	Ra < 0.8 μm (N6)

For more specifications see www.hamiltoncompany.com

Ordering Information

243230 243270		tex C Plas								
	Code									
	1	Material (wetted parts) PVDF / Stainless Steel 2.4602								
	2	PEEK								
	0	special								
		Code Sealing Material (wetted sealings)								
		1		FDA USP		0 ,				
		2	FKM							
		3	FFKM							
			0	special						
			Code	Sensor						
			1	325mm	n PG13,5					
			0	special						
				Code	Process	s Connect	tion			
				1	Flange D	N50				
				2	Flange ANSI 2" special					
				0						
					Code		g Connection			
					1		read female			
					2		ead female			
					0	special				
						Code	Position switch			
						1	pneumatic / without for manual			
						0	special			
Ref							← Order Code			

MasterFit



The MasterFit is a housing for pressurizable pH sensors like the ChemoTrode types. The pressurization ensures a constant outflow of electrolyte. This helps to prevent clogging of the diaphragm and poisoning of the electrolyte. The MasterFit can be used in a huge variety of applications mainly in the chemical industry.

The pressure inside the MasterFit can be controlled via a built-in manometer. Furthermore the liquid level of the electrode can be controlled through the coated glass body of the armature at any time.

Benefits

 $\langle \epsilon_x \rangle$

- ► Sealing feature prevents loss of pressure caused by soiling
- ► Pressure reduction on disassembly
- ► Various o-ring positions available
- ► Easy maintenance

Ordering Information


Туре	Process Connection	Ref
MasterFit 120	G 11/4	237200-OP
MasterFit 150	G 1¼	237225-OP
MasterFit 250	G 1¼	237245-30



For more specifications see www.hamiltoncompany.com

Dimensional drawings / MasterFit 120

all dimensions in mm



Туре	A (armature insertion depth)	B (total length)
MasterFit 120	40 mm	475 mm
MasterFit 150	70 mm	505 mm
MasterFit 250	170 mm	605 mm

Accessories



- Pressure Adapter Ref 237252
- Service Kit for MasterFit Ref 237229
- FFKM Kit for MasterFit Ref 237319

• Flange Adapter for MasterFit* Ref 237910

Safety Socket see page ▶ 152

*The Flange Adapter is used with a MasterFit 120 and a sensor with a shaft length of 150 mm

RetractoMaster



The RetractoMaster is a retractable housing for pressurizable sensors like the ChemoTrode. The pressurization ensures a constant outflow of electrolyte. This helps to prevent clogging of the diaphragm and poisoning of the electrolyte. It allows the operator to mount and dismount the sensors while the process is running. The pressure inside the RetrctoMaster can be controlled via a built-in manometer. Only one press on the red button is needed to move the sensor into or out of the process. Safe sensor handling during process is guaranteed because insertion into a vessel without the sensor is impossible so is removal while in measuring position. O-rings can easily be replaced without special tools. Furthermore the liquid level of the electrode can be controlled through the coated glass body of the armature at any time. The RetractoMaster can be used in a huge variety of applications mainly in the chemical industry.

Benefits

- ► Sensor can be withdrawn from the process for cleaning, calibration or replacement
- ► Easy maintenance
- ▶ Long life time of the sensor due to pressurization of the sensor and the possibility to remove it while the process is running.
- ➤ 3.1 certificate included

Ordering Information

Туре	Process Connection	Ref
RetractoMaster	G 11/4	237250


Specifications Wetted parts Stainless Steel 1.4571 O-ring material FKM O-ring position 22.5 mm Pressure range (relative to ambient) 0 to 6 bar Temperature range -10 to 130 °C Sensor a-length 250 mm Surface finish Ra < 0.4 μm (N5)</td>

For more specifications see www.hamiltoncompany.com

Dimensional drawings

all dimensions in mm

Maintenance position



Measuring position

Accessories



• Pressure Adapter Ref 237252

Safety Socket see page ▶ 152

Safety Socket




The Safety Sockets are hygienic weld-in sockets suitable for hygienic armatures like the FlexiFit Bio. They are available for 3 different o-ring positions to cover different standards. Furthermore you can choose between two kinds of stainless steel and two different angles.

The Safety Socket narrows at the o-ring positions and it seals only if the o-ring of the armature is exactly at the right place. If the process is under pressure, a dripping process medium can be a strong hint that the armature should not be loosened entirely. Therefore the Safety Socket is suited for a wide variety of applications and installations.

Benefits

- ► Safety design, leakage before total release of the armature
- Hygienic surface finish
- ▶ 3 different o-ring positions and two different stainless steels available



SpecificationsWetted partsStainless Steel 1.4435, 1.4404 or 1.4571O-ring material for blind plugEPDMPressure range (relative to ambient)0 to 50 bar (relative to ambient)Temperature range-30 to 160 °CProcess connectionG 1¼Surface finishRa < 0.4 μm (N5)</td>


For more specifications see www.hamiltoncompany.com


Ordering Information

Туре	Steel	Angle	OP	Ref
Safety Socket	1.4404	15	25	242570
Safety Socket	1.4404	15	50	242571
Safety Socket	1.4404	15	55	24257
Safety Socket	1.4404	0	25	24257
Safety Socket	1.4404	0	50	242574
Safety Socket	1.4404	0	55	24257
Safety Socket	1.4435	15	25	24257
Safety Socket	1.4435	15	50	24257
Safety Socket	1.4435	15	55	24257
Safety Socket	1.4435	0	25	24257
Safety Socket	1.4435	0	50	24258
Safety Socket	1.4435	0	55	24258
Weld in socket without safety feat	1.4571 ture	15	25	23720

Dimensional drawings

all dimensions in mm





Accessories



- Blind plug 1.4404-25 Ref 242560
- Blind plug 1.4404-50 Ref 242562
- Blind plug 1.4404-55 Ref 242564
- Blind plug 1.4435-25 Ref 242565
- Blind plug 1.4435-50 Ref 242567
- Blind plug 1.4435-55 Ref 242569
- Blind Plug 1.4571-25 Ref 237230

Hygienic Socket



The Hygienic Socket with its space saving design and simple sterilization is ideal to weld in fermenters or small pipes. The advantages are numerous for many other applications in tanks or pipes for water treatment and in the pharmaceutical and chemical industries.

It is designed for 120 mm sensors and developed for easy installation and maintenance, improve the cleaning process and increase safety. Two "Live Guard" openings provide an indication of sealing failures. The sensor insertion depth can be varied for DO sensors by using the Hamilton DO Adapter.

Benefits


- ► Hygienic design because complete sensor installation with only one wetted o-ring
- Space saving
- ► Cost saving: Socket and Housing all in one
- ► Low maintenance and easy replacement of o-ring




For more specifications see www.hamiltoncompany.com

Dimensional drawings

all dimensions in mm



Hygenic socket with a pH sensor



Hygenic socket with DO adapter and oxygen sensor

Ordering Information

Туре	Ref
Hygienic Socket 1.4404	242535
Hygienic Socket 1.4435	242545
Hygienic Socket 1.4571	242548
Hygienic Socket 2.4602	242550

Accessories



- Hyienic Socket DO Adapter Ref 242538
- Replacment Kit Seal Pusher Ref 242532
- O-ring set EPDM Ref 242595
- **Sensor Dummy 96 mm** Ref 242540
- **Sensor Dummy 117 mm** Ref 242563
- O-ring set FKM Ref 242596
- O-ring set Silicone Ref 242597
- O-ring set FFKM Ref 242598

FlowCell/XL




Hamilton Flow-Through Cells are designed for measuring one or two parameters at a time. Possible combinations are pH/DO and pH/Conductivity. The measurement is done in bypasses when inline measurement is not possible due to small pipe dimensions. Application fields are biotechnology, water treatment and power plants, where reliable measurements have to be carried out in ion-weak media. There are two different sizes of the flow cells available.

Benefits

- ► Flexible design for one or two measuring points
- ▶ PEEK insert of high chemical resistance
- ► Low dead volume
- ➤ Self draining
- ► Internal aseptic clamp pipe connection

Dimensional drawings / 242590

all dimensions in mm



Stainless Steel 1.4435, PEEK
EPDM
0 to 16 bar
-10 to 140 °C
PG 13.5
120 mm
TC 25, TC 50, Swagelok

For more specifications see www.hamiltoncompany.com

Ordering Information

242585						
	Code	Measuring position				
	1	only pH	or Conduc	cell UPW		
	2	only Conductivity or Oxygen pH and Conductivity or Oxygen				
	3					
	4 Conductivity and Oxygen					
	0	special				
		Code	Pipe Co	onnection		
		1	TC25 1/4	"		
		2	TC25 3/8	, ,		
		3	TC25 ½	33		
		4	Swagelo	ok 6 mm		
		5	Swagelo	ok 10 mm * ok 1⁄4"		
		6	Swagelo			
		7	Swagelo	ok 3/8" *		
		8	Swagelo	ok ½" *		
		0	special			
			Code	o-ring material		
			1	EPDM		
			2	FFKM (two measuring positions)		
			3	FFKM (one measuring position)		
	+	+	0	special		
242585 -				← Order Code		

242590							
	Code	Measur	ing position				
	1	only pH	or Conduc	cell UPW			
	2	only Cor	nductivity o	or Oxygen			
	3	pH and	Conductivi	ity or Oxygen			
	4	Conductivity and Oxygen					
	0	special					
		1 2	Pipe Connection TC50 ¾" TC50 1"				
		3	TC50 1.				
		0	special				
			Code	o-ring material			
			1	EPDM			
			2	FFKM (two measuring positions)			
			3	FFKM (one measuring position)			
	*	*	0	special			
242590 -				← Order Code			

*Not self draining



Accessories

- O-ring kit Flow Cell Ref 237387
- O-ring kit Flow Cell XL Ref 237390
- Sensor Dummy 96 mm Ref 242540 117 mm Ref 242563

FlexiFlow SL 10



The FlexiFlow is a flow-through cell. It can be used in all cases where pH or oxygen must be reliably measured in ion-weak media including coolant piping in power generating stations.

The sample is fed into the cell from the bottom at a low flow speed, and out of the cell again at the side. A groove cut into the FlexiFlow allows it to easily be attached anywhere with commercially available screws.

Benefits


- ► Compact design
- Easy to attach to a plate
- ► For use in small pipes where sensors cannot be inserted directly
- ▶ Self draining




For more specifications see www.hamiltoncompany.com

Dimensional drawings

all dimensions in mm





Ordering Information

Ref
237340

SENSOR COMPARISON

Sensor Comparison

pH or ORP sensor

	pH glass type	Nominal measurement range	Recomm. measurement range	Reference system	Reference electrolyte	Diaphragm type	Recomm. min. conductivity (µS/cm)	Nominal temperature range (°C)	Recomm. temperature range (°C)	Nominal pressure max. (bar)	Upside down Installation	Comments
ChemoTrode	PHI	0 to 14	0 to 13	Everef-F	3M KCI-LR	HP ceramic	20	0 to 130	5 to 130	6	No	
ChemoTrode Bridge	PHI	0 to 14	0 to 13	Everef-B	Skylyte	HP ceramic	20	0 to 130	5 to 130	6	No	
ChemoTrode P PHI	PHI	0 to 14	0 to 13	Everef-F	Protelyt	HP ceramic	20	0 to 130	5 to 130	6	No	
FermoTrode	PHI	0 to 14	0 to 13	Everef-F	Skylyte	Coatramic	20	0 to 130	5 to 130	4	No	
EasyControl	HF	0 to 14	0 to 13	Ag/AgCl	Viscous 3M KCI	Ceramic	20	0 to 60	0 to 60	2	No	
InchTrode N100F	HF	0 to 14	2 to 12	Everef-L	Polisolve	Single Pore ring	5	-10 to 130	5 to 100	6	No	
InchTrode N75F	HF	0 to 14	2 to 12	Everef-L	Polisolve	Single Pore ring	5	-10 to 130	5 to 100	6	No	
InchTrode N75FC10	HF	0 to 14	2 to 12	Everef-L	Polisolve	Single Pore ring	5	-10 to 130	5 to 100	6	No	
InchTrode N75P	PHI	0 to 14	2 to 12	Everef-L	Polisolve	Single Pore ring	5	0 to 130	5 to 100	6	No	
InchTrode N75PC10	PHI	0 to 14	2 to 12	Everef-L	Polisolve	Single Pore ring	5	0 to 130	5 to 100	6	No	
IonoTrode	F	0 to 14	0 to 13	Everef	3M KCI	Sleeve	0.2	-10 to 40	-10 to 40	0.5	No	
LIQ-Glass PG	F	1 to 12	1 to 12	Everef	3M KCI-LR	Ceramic	2	-5 to 60	-5 to 60	2	No	
MecoTrode	Н	0 to 14	0 to 14	Everef	Viscous 3M KCI	HP ceramic	50	0 to 130	0 to 130	6	No	0 to 16 bar at 25 °C, 0 to 6 bar at 130 °C
Polilyte Pro	HF	0 to 14	2 to 12	Everef-B	Polisolve	Single Pore	5	-10 to 60	-5 to 60	6	Only VP	
Polyplast Pro	V	0 to 14	2 to 12	Ag/AgCl	Polisolve	Single Pore	50	-10 to 40	0 to 40	6	No	
Polilyte Plus XP	Н	0 to 14	2 to 12	Everef-L	Polisolve Plus	Single Pore	5	0 to 130	0 to 130	16	Only VP	0 to 50 bar (60 °C), 0 to 20 bar (100 °C), 0 to 16 bar (130 °C)
pH families												
Polilyte Plus H	Н	0 to 14	2 to 12	Everef-L	Polisolve Plus	Single Pore	5	0 to 130	0 to 130	10	Only VP/MS	Predecessor: Polilyte Plus, Polilyte HT
Polilyte Plus HB	НВ	0 to 14	2 to 12	Everef-L	Polisolve Plus	Single Pore	5	0 to 130	0 to 130	10	Only VP/MS	
Polilyte Plus HF	HF	0 to 14	2 to 12	Everef-L	Polisolve Plus	Single Pore	5	-10 to 100	-10 to 100	16	Only VP/MS	Predecessor: ClaryTrode
Polilyte Plus PHI	PHI	0 to 14	2 to 12	Everef-L	Polisolve Plus	Single Pore	5	0 to 130	5 to 130	10	Only VP/MS	Predecessor: Polyclave
EasyFerm Plus PHI	PHI	0 to 14	2 to 12	Everef-F	Phermlyte	HP Coatramic	100	0 to 140	5 to 140	6	No	
EasyFerm Plus HB	НВ	0 to 14	2 to 12	Everef-F	Phermlyte	HP Coatramic	100	0 to 140	5 to 140	6	No	
EasyFerm Bio PHI	PHI	0 to 14	2 to 13	Everef-F	Foodlyte	HP Coatramic	100	0 to 140	0 to 140	6	No	
EasyFerm Bio HB	НВ	0 to 14	2 to 13	Everef-F	Foodlyte	HP Coatramic	100	0 to 140	0 to 140	6	No	
ORP Sensors												
ChemoTrode ORP	Platinum ring	± 2000 mV	± 2000 mV	Everef-F	3M KCI-LR	HP ceramic	20	0 to 130	0 to 130	6	No	
EasyControl ORP	Platinum wire	± 2000 mV	± 2000 mV	Ag/AgCI	Gel	Ceramic	20	0 to 60	0 to 60	2	No	
OxyTrode Pt	Platinum wire	± 2000 mV	± 2000 mV	Everef	Viscous 3M KCI	HP ceramic	50	0 to 130	0 to 130	6	No	
Polilyte RX	Platinum wire	± 2000 mV	± 2000 mV	Everef-B	Polisolve	Single Pore	5	-10 to 60	-10 to 60	6	No	
Polyplast Pro RX	Platinum wire	± 2000 mV	± 2000 mV	Ag/AgCI	Polisolve	Single Pore	50	-10 to 40	-10 to 40	6	No	
EasyFerm Plus ORP	Platinum wire	± 2000 mV	± 2000 mV	Everef-F	Phermlyte	HP Coatramic	100	0 to 140	5 to 140	6	No	Arc: ± 1500 mV
Polilyte Plus ORP	Platinum ring	± 2000 mV	± 2000 mV	Everef-L	Polisolve Plus	Single Pore	5	0 to 130	0 to 130	10	Only VP	Arc: ± 1500 mV, 0 to 16 bar at 100 °C, 0 to 3 bar at 140 °C

SENSOR COMPARISON
SENSOR COMPARISON

DO sensor

	Measurement principle	Nominal measurement range (DO)	Nominal temperature range	Measurement temperature range	Nominal pressure max. (bar)	Compatible caps / membranes
VisiFerm RS485	Optical	4 ppb to 25 ppm	-10 to 140 °C	-10 to 85 °C	12	H0, H2, H3, H4
VisiFerm mA	Optical	4 ppb to 25 ppm	-10 to 140 °C	-10 to 85 °C	12	H3, H4
VisiTrace mA	Optical	1 ppb to 2 ppm	-10 to 140 °C	-10 to 85 °C	12	L1
VisiWater DO P Arc	Optical	0 to 40 ppm	0 to 60 °C	0 to 60 °C	12	H20
OxyFerm FDA	Amperometric	10 ppb to 40 ppm	0 to 130 °C	0 to 130 °C	4	FDA, CIP, standard
OxyGold B	Amperometric	8 ppb to 40 ppm	0 to 100 °C	0 to 100 °C	12	OxyGold
OxyGold G	Amperometric	1 ppb to 40 ppm	0 to 130 °C	0 to 130 °C	12	OxyGold
Oxysens	Amperometric	40 ppb to 40 ppm	0 to 60 °C	0 to 60 °C	4	none

Conductivity sensor

Conducell 4US 4 pole contacting 0.1 μS/cm to 500 mS/cm -20 to 135 °C 0.147/cm 6 Stainless steel 1.4435 Conducell UPW 2 pole contacting 0.01 to 1500 μS/cm 0 to 130 °C < 0.1/cm 10 Stainless steel 1.4435							
Conducell 4US 4 pole contacting 0.1 μS/cm to 500 mS/cm -20 to 135 °C 0.147/cm 6 Stainless steel 1.4435 Conducell UPW 2 pole contacting 0.01 to 1500 μS/cm 0 to 130 °C < 0.1/cm 10 Stainless steel 1.4435		Measurement principle	Nominal measurement range	Nominal temperature range	Cell constant	Nominal pressure max. (bar)	Electrodes materials available
ConduceII UPW 2 pole contacting 0.01 to 1500 μS/cm 0 to 130 °C < 0.1/cm 10 Stainless steel 1.4435	Conducell 4UxF	4 pole contacting	1 μS/cm to 300 mS/cm	-20 to 150 °C	0.36/cm	20 (135 °C)	Stainless steel 1.4435, Titanium, Hastelloy C 2.4602, Platinum
	Conducell 4US	4 pole contacting	0.1 μS/cm to 500 mS/cm	-20 to 135 °C	0.147/cm	6	Stainless steel 1.4435
Conducell 2DC-PG 2 pole contacting 10 uS/cm to 20 mS/cm -5 to 80 °C 1/cm 6 Graphite	Conducell UPW	2 pole contacting	0.01 to 1500 μS/cm	0 to 130 °C	< 0.1/cm	10	Stainless steel 1.4435
	Conducell 2DC-PG	2 pole contacting	10 μ S/cm to 20 mS/cm	-5 to 80 °C	1/cm	6	Graphite

Safety First

Hamilton Offers More Certificates Then Ever


Many industrial processes are in hazardous environments and require suitable equipment with the European ATEX, the British UKEX or the global IECEx approval. Hamilton provides safe sensors and housings since many years for these applications. In case a gas atmosphere and a dust atmosphere are or could be present at the same time, the risk of explosion must be examined carefully and special precautions may be necessary. Typical gas atmospheres can be found in oil refineries, printing industries and biogas plants. Dust atmospheres can be found in underground coalmines, woodworking areas and in all kind of mills. In the chemical industry both atmospheres can be found.

ATEX is the widely used synonym for the ATEX directives of the European Union. ATEX stands for the French abbreviation «ATmosphère EXplosible». The objective of ATEX is to ensure the free movement of goods throughout the European Union, by offering one harmonized compliance procedure accepted by all EU countries. This means that different national standards within the EU are obsolete. ATEX covers equipment only. Equipment for hazardous areas requires an ATEX approval when sold within the European Union.

The **UKEX** regulation applies to Great Britain and corresponds to the ATEX directive.

The **IECEx** system is a conformity assessment system of the International Electrical Commission (IEC). It is the objective of the IECEx system to facilitate international trade in equipment and services. Currently Australia, New Zealand, and Singapore accept the IECEx certificate of conformity as meeting all of the national requirements for Ex Certification. No further national certification is required. The IECEx is also accepted in many other countries.

Marking sensors or housings for ATEX / IECEx is as follows:



Example OxyFerm FDA

The temperature value x in dust atmospheres needs to be calculated.

The table gives an overview of the approvals available for the different product lines. Detailed information about a specific product can be found on the Hamilton website their spec sheets. For general overview please refer to:

www.hamiltoncompany.com/support/process-analytics/certificates/products-for-explosive-atmospheres

	ATE		UK	EX	IE	CEx
Sensor/Housing	Gas	Dust	Gas	Dust	Gas	Dust
Analog Sensors	/	/	/	/	/	/
Housings	/	/	/	/	/	/
Arc	_	_	-	-	-	_
Memosens	/	_	/	_	/	_
VisiFerm mA	/	/	/	/	/	/
VisiTrace mA	/	/	/	/	/	/

Alphabetical Index

2-pole conductivity sensor	62	FermoTrode	30, 160	ORP buffers	101	Retractex C Plastic	142, 14
3 M KCI	105	FFKM Kit MasterFit	149	OxyFerm	92	Retractex C Plastic (pneuma	atic) 14
4-pole conductivity sensor	60	FFKM Kit RetractoFit	131	OxyFerm CIP	93	Retractex C Plastic LT (pneu	matic) 14
,		Fix cable	35, 61, 91, 117	OxyFerm FDA	92, 93, 104, 108, 112	Retractex C Plastic LT M (ma	anual) 14
A Arc	12, 13	Flange Adapter MasterFit	149	OxyFerm Membrane Kit	104	Retractex C Plastic M (manu	
Arc Accessories	117	FlexiFit	128	Oxygen Accessories	104	Retractex C Steel	140, 14
Arc ECS Adapter	15, 115	FlexiFit Bio	128	OxyGold B	94, 162	Retractex C Steel (pneumati	c) 14
Arc Module Cond-P SU	59	FlexiFit TC	128	OxyGold G	96, 162	Retractex C Steel LT	14
Arc Module Cond-1 30	71	FlexiFit U Bio	128	OxyGold Membrane Kit	104	Retractex C Steel LT (pneun	
		FlexiFit VV-0		Oxydoid Membrane Kit Oxylyte	104, 105	Retractex C Steel LT (prieding	
Arc Module pH SU	25		128				
Arc Modbus OPC Converter	118	FlexiFit VV-15	128	Oxylyte B	105	Retractex C Steel M (manua	
Arc Sensor USB Power Cable	116	FlexiFlow SL 10	158	Oxylyte G	105	RetractoFit	13
Arc Systems	15	Flow Cell	156	Oxysens	98, 162	RetractoFit Bio	13
Arc View Mobile	117	Foodlyte	8, 22, 23, 160	OxyTrode Pt	48, 160	RetractoFit Bio 25	13
Arc Wi 1G Adapter BT	117					RetractoFit Bio 55	13
Arc Wi 2G Adapter BT	117	H H100 Cond	124	P Panel-mount Kit	121, 122	RetractoFit PEEK 25	13
Autoclavation Cap OxyFerm	104	H100 Condl	122	рН	16-39	RetractoMaster	15
		H100 DO	123	pH buffers	10, 100		
B Blind plug	153, 156	H100 pH	120	pH glasses	8	S S7	29, 31, 33, 47, 10
Bluetooth	88	H220X	124	pH pathfinder	17	S8 19, 21, 2	23, 25, 37, 39, 43, 45, 49, 51, 53, 10
BNC plug	110, 115	Hamilton Customized Products	119	pH Port	27	Safety Socket	15
2.10 p.a.g		Housings	126-159	Phermlyte	20, 21, 44, 160	Sensor Cable M12	10
C Cables	108	HyCIP	134	Pipe-mount Kit	120, 121, 122, 123	Sensor Power Cable M12	114, 11
Calibration station	77, 87, 89	Hygienic Socket	154	Polarization Module	104	Service Kit FFKM	12
Cell Density	4, 66-73	Trygleriic Socket	154	Polarization Module B	104	Service Kit FKM	12
ChemoTrode	28, 46, 160	I InchTrode	34, 160	Polarization Module G	104	Service Kit FlexiFit Bio	12
				Polarization Module T	104	Service Kit Hexil it Bio	14
ChemoTrode Bridge	28, 160	Incyte Arc	4, 70	Polarization Module 1 Polityte Plus		Service Kit MasterFit Service Kit RetractoFit	13
ChemoTrode ORP	46	Incyte SU	5, 72		18, 42, 160		
ChemoTrode P	28	Insertion tube	136, 138, 140, 142, 144, 146	Polilyte Plus ORP	42, 160	Single Pore	9, 19, 35, 39, 43, 51, 16
Clark cell	79	IonoTrode	32, 160	Polilyte Plus XP	18, 19, 160	Skylyte	29, 31, 105, 16
Cleaning solution set	105			Polilyte Pro	36, 160	Storage solution	10
Conducell 2DC-PG	64, 162	K K8	21, 23, 27, 108, 111	Polilyte RX	50, 160	System Installations	14, 1
Conducell 4US	60, 162			Polisolve	35, 36, 37, 38		
Conducell 4USF	56, 162	L Lemo plug	110, 113	Polisolve Plus	9, 18, 19	T T82/D4	97, 113, 11
Conducell 4UxF	56, 111, 162	Liq-Glass PG	38, 160	Polyplast Pro	36, 160	Transmitter	120-12
Conducell SU	58	Liquid Earth	20, 22	Polyplast Pro RX	50, 160	Triclamp (TC)	57, 61, 63, 67, 128, 135, 141, 15
Conducell UPW	62, 113, 162			Pressure Adapter	149, 151		
Conductivity	55	M M12	108, 113, 114, 116	Protective Hood	120, 121, 122, 123	V Varivent (VV)	56, 128, 13
Conductivity Standards	10, 102	MasterFit	148	Protelyte	29, 105	VisiFerm DO	1
CO ₂ NTROL	6, 74-77	MecoTrode	24, 160	PTFE scraper	136, 138, 140, 142, 144, 146	VisiFerm DO SU	5, 8
Customized Products	118	Membrane Kit CIP	104			VisiFerm mA	84, 16
		Membrane Kit FDA	104	R Replacement Cathode OxyFerm	93, 104	VisiFerm RS485	6, 80, 16
D DCO ₂	6. 74-77	Memosens (MS)	19, 21, 23, 106, 108, 112, 124	Replacement Cathode OxyGold B	95, 104	VisiTrace mA	88, 16
Dencytee Arc	4, 72	Modbus Profibus Converter	117	Replacement Cathode OxyGold G		VisiTrace RS485	6, 86, 16
	78-99	Modbas i Tolibas Colivertei	117	Replacement Kit Seal Pusher	155	VisiWater DO P	90, 16
DO, Dissolved Oxygen		ODO Cap H0 / H2	01	Retractex	134-147	VP	108, 10
DO Adapter	154, 155		81	Retractex B		VP 6	11
DuraCal pH buffers	10, 100	ODO Cap H3 / H4	85		134		
		ODO Cap S0 / S2	83	Retractex B (pneumatic)	135	VP 8	11
E EasyControl	38, 52, 160	OneFerm pH	5, 26	Retractex B M (manual)	135		
EasyControl ORP	52, 160	Optical oxygen sensor	80-89	Retractex BC	136-139	W Weld-in socket	15
EasyFerm Bio	22, 160	O-ring Kit Flow Cell	157	Retractex BC Plastic	138	Wireless Converter BT	11
EasyFerm Plus	20, 44, 160	O-ring Kit Flow Cell XL	157	Retractex BC Plastic (pneumatic)	139		
EasyFerm Plus ORP	44, 160	O-ring set EPDM	155	Retractex BC Plastic M (manual)	139		
Electrochemical oxygen sensor	92-99	O-ring set FFKM	155	Retractex BC Steel	136		
Electrolyte	104, 105	O-ring set FKM	155	Retractex BC Steel (pneumatic)	137		
<u> </u>		O-ring set Silicone	155	Retractex BC Steel M (manual)	137		
F FDA Service Kit	133	ORP	40-53	Retractex C	140-147		
	100	·	10 00	-	· · · · · · · ·		

BioConnect® is a registered trademark of NEUMO GmbH + Co. KG, Knittlingen (D).

Bluetooth® is a registered trademark of Bluetooth SIG Inc., Kirkland WA, USA.

Knick® is a registered trademark by KNICK Elektronische Messgeräte, Berlin (D).

Memosens® is a registered trademark of Endress + Hauser, Reinach (D).

Tuchenhagen Varivent® is a registered trademark of GEA Tuchenhagen GmbH.

Unigate® is a registered trademark of Deutschmann Automation GmbH & Co. KG, Bad Camberg (D).

App Store, iOS, and their logos are registered trademarks of Apple Inc. in the US and other countries.

Android, Google Play, and their logos are registered trademarks of Google Inc. in the US and/or other countries.

Windows and their logos are registered trademarks of Microsoft Corporation in the US and other countries.

All other trademarks are owned and/or registered by Hamilton Bonaduz AG.

 $\ensuremath{\texttt{©}}$ 2022 Hamilton Bonaduz AG. All rights reserved.

REF 690127/13 - $^{\text{M}}$ 07/2022



HAMILT®N®

Web: www.hamiltoncompany.com USA: 800-648-5950 Europe: +41-58-610-10-10

Hamilton Americas & Pacific Rim Hamilton Company Inc. 4970 Energy Way Reno, Nevada 89502 USA Tel: +1-775-858-3000 Fax: +1-775-856-7259 sales@hamiltoncompany.com

Hamilton Europe, Asia & Africa Hamilton Bonaduz AG Via Crusch 8 CH-7402 Bonaduz, Switzerland Tel: +41-58-610-10-10 contact.pa.ch@hamilton.ch

To find a representative in your area, please visit www.hamiltoncompany.com.