

CAT-01 CAT-03

SPECIFICATION SHEET

Product : Cation filter
Model : Cat-01 Cat-03

Filter length : 630 mm *
Filter contents : 1.0 / 3.0 Ltr *

Max. Sample pressure : 3,5 Bar abs (1 liter version Cat-01)
Max. Sample pressure : 4,0 Bar abs (1 liter version Cat-03)
Max. Sample temp. : 70 °C
Max. Sample flow : 0,5 ltr/minute (1 liter version Cat-01)
Max. Sample flow : 1,2 ltr/minute (3 liter version Cat-03)

Typical applications: Boiler IP and HP Steam multiple streams, Process cooling of liquids

Features:

- Rugged design for an industrial environment.
- 100% chemical resistant, no discoloration of tube.
- Connection allows easy access for filter inspection.
- Multi Instruments Analytical B.V. always uses regenerative exchange resins with a colour indicator.
- The Cation Filters exchange vessels can be easily changed and facilitate fast and easy maintenance of the system.
- Compact body size.

* Other contents/dimensions on request

