


iFiD Master

Remote Control

Control via TCP/IP




Product description

The iFiD Master is an operating software for your FID. It is controlled via a TCP/IP network connection, which is possible from any Windows-compatible device.

You can easily operate the FID from the office or operate the FID completely from the office or from home.

Special advantages

- Control via TCP/IP
- Data recording directly to the PC, tablet or laptop
- Full access to your FID from the office or from home
- Clear and user-friendly user interface
- Quick reaction and response from the FID
- Reliable connection
- Remote diagnosis for problems
- Extended additional function of the operation
- Integrated query commands



Operating principle

Connect the FID to a network and connect it to the PC, tablet or laptop via the shared port and IP address.

